Nuclear Magnetic Resonance Facility

User's Manual

Department of Chemistry and Biochemistry University of Oklahoma

Version 2.1.03 March 16, 2011 Stephenson Life Center, Room 1700

Acknowledgements

This manual was constructed using VNMRJ2.2C and VNMRJ2.2D with a Red Hat Linux 5.1 operating system on the computers within the NMR facility at the University of Oklahoma. Material was also incorporated from the Structural Biology NMR lab User's Manual from the University of Minnesota with the express permission and review of Dr. Beverly Ostrowski.

The NMR Facility of the Department of Chemistry and Biochemistry at the University of Oklahoma would like to express gratitude to the Structural Biology NMR Resource at the University of Minnesota for their willingness share their knowledge, resources and training.

Suggested Reference Materials

Books

Modern NMR Spectroscopy by Sanders and Hunter High-resolution NMR techniques in organic chemistry by Claridge, Timothy D. W Spectrometric Identification of Organic Compounds by Silverstein and Webster 200 and More NMR Experiments by Berger, and Braum

Article

Reynolds, W.F.; Enriquez, Raul G.; "Choosing the Best Pulse Sequences, Acquisition Parameters, Postacquisisiton Processing Strategies, and Probes for Natural Product Structure Elucidation by NMR Spectroscopy", J. Nat. Proc., **2002**, 65, 221-224.

Web Resource Basic NMR Concepts: A Guide for the Modern Laboratory by Dr. Daniel Holmes. http://www.chemistry.msu.edu/facilities/nmr/

	Contents							
	Abbreviations and Terms	5						
	NMR Laboratory Policies	6						
	NMR Training Policies	7						
Ι.	Preparing Samples	7						
	Tubes	7						
	Handling and Preparation	8						
	Sources of Contamination	8						
П.	Equipment	9						
	Spectrometers	9						
	Probes	10						
III.	Choosing the Spectrometer	10						
	Instrument specifications and set-up	10						
	Sensitivity	11						
	Resolution	12						
	Dynamic Molecules	12						
IV.	Software	13						
	Intro to VNMRJ	13						
V.	Data Files	15						
	Creating Data Directories	15						
	Creating a Desktop Shortcut to the Data Directory	16						
	Creating Data Transfer Shortcuts on the Red Hat Desktop	17						
	Transferring the data	19						
VI.	Probe Files	21						
	Instrument Calibration	21						
	Probe Files	21						
	Changing the Probe File	22						
VII.	Gradient Shimmaps	23						
VIII.	Using the Varian Mercury VX-300 MHz Spectrometer	25						
	300 NMR sign-up rules	25						
	Mercury 300 Tests and Assignment Certification	26						
	Collecting a 1D proton on the Mercury VX-300 MHz NMR	27						
IX.	Using the Varian VNMRS-400 MHz Spectrometer	33						
	400 NMR sign-up rules	33						
	VNMRS- 400 Tests and Assignment Certification	34						
	Collecting a 1D proton on the VNMRS-400 MHz NMR	35						
	Tuning the broadband probe using protune('calibrate')	42						
Χ.	Using the Varian VNMRS-500 MHz Spectrometer	44						
	500 NMR sign-up rules	44						
	VNMRS- 500 Tests and Assignment Certification	45						
	Tuning the 500 probes	46						
	Using mtune	47						
XI.	Changing Basic 1D parameters	48						
	Setting the frequency window (adjusting sw and tof)	49						
	Setting the first delay (d1) and acquisition time (at)	50						

XII.	Basic 1D processing	51
	Opening Saved Data	51
	Toolbar Icons	53
	Phasing	54
	Referencing the Spectrum	55
	Weighting Functions	56
	Zero Filling	58
	Measuring Signal-to-noise	58
	Finding digital resolution	59
	Integration	59
	Baseline Correction	62
	Peak Picking	63
	Basic 1D Plotting	64
	Making PDF Files	64
	Converting 1D Data to ascii Format	65
XIII.	Finding a 90 degree pulse width	66
XIV.	Basic 13C 1-D NMR	71
XV.	2D- Experiments	74
	2D Parameter Sets	74
	General Considerations for doing 2D Experiments	75
	COSY	76
	NOESY	78
	HSQC	80
	HMBC	81
XVI.	2D Processing and Printing	83
	2D Tool Bar	83
	Manipulating 2D Data	83
	Processing 2D data sets manually	84
	Manual Phase correction	85
	Referencing the 2D data	86
	Plotting the 2D data	87
	Appendix	88
A.1	Mercury VX-300 Practical Assignment #1	88
A.2	Mercury VX-300 Practical Assignment #2	89
A.3	VNMRS-400 Practical Assignment #1	91
A.4	VNMRS-400 Practical Assignment #2	96
A.5	VNMRS-500 Equipment components	99
A.6	Transferring Data to your pc using Winscp	105
A.7	Remotely accessing the 400 NMR via PUTTY and REALVNC	107
A.8	Glossary of Common NMR Commands and Terms	113
A.9	A few useful Linux commands.	120

Abbreviations and Terms

Correlation Spectroscopy
Double-quantum filtered correlation spectroscopy
"Indirectly" detected dimension
"Directly" detected dimension
Free Induction Decay
Heteronuclear Multiple Bond Coherence
Heteronuclear Multiple Quantum Coherence
Heteronuclear Single Quantum Coherence
Hertz
Megahertz
Molecular Weight
Nuclear Magnetic Resonance
Nuclear Overhauser Effect
Nuclear Overhauser Effect Spectroscopy
Pulsed field gradient
parts per million
part of the spectrometer where RF signal is sent and received
used interchangeably for 'H
radio frequency
Rotating frame Overhauser effect spectroscopy
Longitudinal Relaxation
Transverse Relaxation
Total Correlation Spectroscopy
Varian's software used to interface with console
variable temperature control
used to refer to anything not 'H

NMR Laboratory Policies

Laboratory Information NMR Facility Manager: Office: Phone: Email:

Dr. Susan L. Nimmo SLSC 1700 627-7044 SusanA@ou.edu

Instrument Problems

Please contact Dr. Nimmo by cell phone and email with details of the problem.

Instrument Log Book

Please record your name, group, general experimental details and any problems that are encountered when you use the instrument.

Sample Breakage

In the event of a sample tube breaking in the magnet please contact the NMR facility manager immediately. This should be done both by phone and email. It is imperative that no one use the instrument until after Dr. Nimmo has accessed the situation. The student should place a "broken sample" note on the keyboard and lock the computer screen.

Faces On-Line Signup

The faces on-line scheduling program hosted by the Complex Carbohydrate Research Center at the University of Georgia is used for instrument sign-up for all instruments our laboratory. The page is found at <u>http://faces.ccrc.uga.edu/</u> It can be accessed from any computer which is attached to the internet. Individual login and passwords will be set up for users at the initial training sessions. You will only be allowed to sign-up for instruments that you have been trained to use. Please read the news page of the sign-up for instrument status and special instructions for each instruments. All users will be held responsible for following specific instructions found on this news page.

Funding Acknowledgements

Please acknowledge the instrument funding and notify the NMR facility manager by email when publishing data collected on the 300 and 400 spectrometers. It is important for us to keep statistics for both past and future proposals. The Varian Mercury VX-300 NMR Spectrometer was purchased in 2000 by the multi-user NSF grant CHE#0077707. The Varian VNMRS-400 NMR Spectrometer was purchased in 2007 by the multi-user NSF grant CHE#0639199

NMR Training Policies

All students will begin training on the 300 NMR. Training sessions will be approximately 45 minutes in duration and will be held as many times as necessary for the student to gain competency. More advanced students and post-docs may gain certification by demonstrating their ability to use the instruments and by completing the assignments.

Certification to use the 300 NMR will be given as NMR Test #300-1, Test #300-2 and Practical Assignments #300-1 and #300-2 are completed as approved by Dr. Nimmo.

Training will begin on the 400 upon successful completion of 300 training. Certification to use the 400 NMR will be given as NMR Test #400-1, and Practical Assignments #400-1 and #400-2 are completed as approved by Dr. Nimmo.

Training will begin on the 500 upon successful completion of 400 training and Test #500-1. Certification to use the 500 NMR between 8-5 Monday-Friday will be given as NMR Test #500-2 is completed as approved by Dr. Nimmo. Unlimited access will be given upon the completion of NMR Test #500-3

Additional training may be requested on an individual basis. But please understand that training sessions *will not* last longer than 1 hour per session.

Variable temperature training will be given as needed. Supervision of variable temperature experiments will continue until both the NMR facility manager and the student are completely confident of the student's ability to safely carry out these experiments independently.

I. Preparing Samples

Tubes

For a standard samples, use the equivalent of 5 mm tubes Wilmad 528-pp-8 (available from the stockroom) (rated 500 MHz and 8 inches long). Eight inch tubes are required for the robotic insertion on the 400. The sample volume should be at least 500 microliters, typically 650-800 microliters. Minimum sample height rule of thumb: coil + 3X diameter of the sample. The probe coils are 16 mm (18mm in the 500 MHz triple resonance and 500 MHz indirect detection probes). Standard NMR tubes are 5mm. Samples shorter than 500 microliters are very difficult to shim. Users running samples containing hazardous materials are encouraged to use thick-walled tubes for safety and should notify facility staff. **Do not use chipped, cracked or scratched tubes**. Labels for NMR tubes can also be purchased from most tube vendors. Alternatives for smaller volume

samples include Shigemi tubes or susceptibility plugs that can be purchased from Wilmad and generally work better for small volumes than thick-walled tubes. Shigemi tubes, 5mm, use volumes of 270-300 microliters and are made of glass that is matched to the magnetic susceptibility of the solvent to be used. Using more than 800 microliters in a regular tube is not recommended because the sample height will extend past the controlled variable temperature range and a temperature gradient can be created in the sample. Some vendors also sell thick-walled tubes that are useful when you have a hazardous sample and are concerned about tube breakage.

Handling and Preparation

Samples should be clear of precipitate and particulates if possible. Many options are available for filtering samples. Centrifugation can also eliminate precipitate from samples. Be sure to wipe off the NMR tube with a Kimwipe around the bottom where the sample may have been touched. Samples require deuterium solvent for a lock signal. Higher quality deuterium solvents, purchased in ampules, will have a smaller residual solvent peak. Samples can be run without deuterium, no lock, for very short runs if necessary. Store deuterated solvents in well-sealed containers on dessicant. Samples should be transferred into and out of NMR tubes using long pipettes that reach the bottom of the NMR tube. Extra long glass pipettes are available from Wilmad. Transfer samples smoothly to avoid losing liquid to the sides of the tube or adding air bubbles to the sample. Water (H_2O) samples require approximately 10% D_2O for a lock signal. Controlling pH in biological samples and in many water samples is highly recommended. pH can change chemical shifts drastically in many aqueous samples. For example, phosphorous shifts are changed dramatically with pH. Cleaning tubes can be accomplished using a tube cleaner which attaches to an aspirator. These are available from several companies. Rinsing the tube several times with water or another solvent if necessary, then with isopropanol or acetone. Do not dry tubes in an oven. The heat in the oven will warp the glass. Blowing nitrogen or other dry gas through the tube will help remove residual solvent, alternatively, pull a vacuum on the tube for several minutes. Another source of cleaning information can be found on the Wilmad website at: http://www.wilmad-labglass.com/services/NMR 010.isp

Sources of Contamination

Tables of chemical shifts of common contaminants are readily available online or in books. Membranes and membranes in concentrators, dialysis bags, buffers (i.e. Tris), detergents, and dirty tubes can all be sources of contamination. Contamination peaks are usually easy to spot because they are often of a different linewidth and different ratio to the sample being observed. Buffers such as Tris contain proton signals that contaminate spectra. Tris is difficult to remove from a protein sample and should be avoided or deuterated if possible. Many commercial membranes are stored in glycerol and should be washed thoroughly before use with a protein sample. Samples prepped in fully deuterated solvent can pick up water from the air very easily and should be sealed tightly. Often the cap can be sealed with a very small piece of Parafilm and samples stored on dessicant in the refrigerator. Samples prepared in "100%" deuterated solvent should be prepared with minimal air exposure or prepared in a dry box or dry bag. Keep bottled solvents on dessicant or, preferably, use sealed ampules of solvent. Solvents from a general use bottle can also become easily contaminated. Occasionally, even newly purchased solvents can contain contamination that is detected when working at low concentrations in high field magnets. To determine if the solvent is the source, run a sample of solvent only. Often a peak near zero is also often seen in samples that can be due to grease used on glassware. Drying tubes with cans of compressed air can also generate a contaminate peak from the propellants and should be avoided. One useful reference in identifying impurities is *J Org Chem*, 62, 7512-7515 (1997), "NMR Chemical Shifts of Common Laboratory Solvents as Trace Impurities".

II. Equipment

Spectrometers

There are three spectrometers located in the NMR laboratory: Varian Mercury VX-300 MHz, Varian VNMRS-400 MHz and VNMRS-500 MHz. All spectrometers are equipped with wave form generators and pulse field gradients. The frequency at which protons precess in the particular strength of magnetic field is used to designate the magnets. For example, protons precess at 300 MHz in the 300 MHz magnet, but the magnet strength is actually 7 Tesla. Further, carbon precesses at about 75 MHz in a 300 MHz magnet. Nuclei such as carbon, nitrogen and phosphorous resonate at much lower frequencies than proton. Often, nuclei that resonant at the higher frequencies such as proton and fluorine are called "high band" and other nuclei like carbon and nitrogen are called "low band". The basic components of the spectrometer include: workstation, console, magnet and probe.

Workstation

The computer workstation is where most of the operation of the instrument occurs including data collection and simple processing. Extensive processing should be done on offline workstations that do not absorb instrument time. The workstation communicates with the console that, in turn, controls the console and the probe in the magnet. The workstations are Dell PC's running under Red Hat Linux 5.1.

Console

The console contains the radio frequency generators, amplifiers, a variable temperature controller, pulsed-field gradient generator, waveform generators, and other computer components. In typical operation, a user will very rarely, if ever, need to interact with the console.

Magnet

The magnetic field in all of the instruments is generated by a current flowing through a solenoid of superconducting wire. For the wire to be superconducting, the wire must stay at liquid helium temperature (4 K) or below. Therefore, the cryostats are filled with liquid helium and outer liquid nitrogen to keep the magnets cold. If the magnets warm up above that temperature, a quench can occur. A quench is when the current in the magnet coil is lost. If a quench occurs, it is usually accompanied by a loud noise followed by fast release of helium gas from the cryostat. If this occurs, please leave the lab as quickly as possible. The magnet is contained inside the silver dewar. The magnets are mounted on vibration legs. The air legs maintain level and stable against small vibrations by air pressure. Therefore, do not lean against the magnets because the magnet will rock. Try to avoid walking around near the magnet during an experiment because it can contribute to vibrations. Typically, the only time a user needs to go near the magnets is to insert the sample and tune the probe. Never take metal or magnetic objects near the magnets. Always check pockets and person for these things before approaching the magnet. Non-digital watches, cards with magnetic strips, and magnetic media (such as disks) will also be affected by the magnetic field. Also inside the magnet are the shim coils. Shim coils are a collection of electrical coils used to remove residual magnet field inhomogeneities. The temperature in the bore of the magnet, where the sample will sit, is controlled by a variable temperature controller and is typically set at 25 deg. C.

Probes

The probe is inside the bore of the magnet. The probe contains the transmitter/receiver coils on where pulses go into the sample and RF frequencies come out. Probes can be changed in a few minutes by facility personnel and have different configurations depending on the application intended. Probes are only changed by facility personnel or by specially trained users. The type of probe selected is determined by the nucleus to be detected and the specific experiment. The 400 and 500 spectrometers have indirect detection probes. These probes are the best choice for direct proton detection or indirect detect experiments. This probe has the proton transmitter/receiver coil closest to the sample and is, therefore, most sensitive for proton detection. The 300 and 400 spectrometers have broadband probes which typically have the X-nucleus coil closer and are, therefore, more sensitive for nuclei like carbon. Carbon can still be detected directly using an indirect detection probe, but it will have a much lower signal-to-noise. Broadband probes or probes that include a X-nucleus can be tuned to a different nucleus depending on the tuning range of the probe. Direct detection of an X-nucleus is best done with a broadband probe.

III. Choosing the Spectrometer

Instrument Specifications and Set-up

It is important to choose the appropriate instrument in the laboratory to answer the relevant experimental questions. Important things to consider include type of probe on the magnet, sensitivity of the system in regards to sample concentration, variable temperature set up, and magnetic field strength.

Varian Mercury 300 MHz NMR Spectrometer

The 300 NMR is especially easy to use. The afternoons are set up in 15 minute time slots because it is possible for anyone to collect a ¹H NMR spectrum in a very short amount of time. The spectrometer has a 4-nuclei auto-switchable probe. This means that ¹H, ¹⁹F, ¹³C, and ³¹P experiments can be collected without changing cables or tuning the probe. This instrument is also capable of running more advanced 1D and 2D experiments, but the sensitivity is less than on the other spectrometers and therefore requires higher sample concentration. The spectrometer is running VNMRJ 2.2D software with chempack 4.1.

Varian VNMRS 400 MHz NMR Spectrometer

The 400 NMR is equipped with a 50 sample autosampler. Therefore this is the best spectrometer for running several samples. It can be set up to automatically run multiple experiments on different samples. It is more sensitive than the 300 NMR and requires less sample concentration. It is equipped with dual broad band probe and automatic tuning. As a result this is the best spectrometer for collecting direct detect broad band spectra. This spectrometer has extended variable temperature accessories and can conveniently collect variable temperature experiments. It also has an indirect detect 1H probe for running indirect detection experiments with lower sample concentration. The spectrometer is running VNMRJ 2.2C software with chempack 4.1.

Varian VNMRS 500 MHz NMR Spectrometer

The 500 NMR is three channel system with both a triple resonance probe and an two channel indirect detection probe. This system has the highest resolution, best sensitivity in the laboratory. Longer 2D and 3D experiments have preference on this instrument over short 1D ¹H experiments. The spectrometer is running VNMRJ 2.2 c software with an chempack 4.1 and Biopack options.

Sensitivity

Signal to noise increases as field strength is increased. Signal to noise is also very dependent on the probe. Indirect detection probes are constructed to maximize the proton sensitivity, while direct detect probes are constructed to maximize broad band signal intensity. The table below lists the signal to noise values for the NMR instruments according to the field strength and probe. These measurements were made after reinstallation in our new laboratory with the Varian standard sensitivity standards.

	1H	F19	13C	P31	N15
300 MHz NMR					
4-nuclei probe	125:1	85:1	90:1	72:1	
400 MHz NMR					
Broad-band probe	300:1	325:1	180:1	180:1	20:1
Indirect Detection Probe	725:1				

500 MHz NMR			
Indirect Detection Probe	965:1	 75:1	
Triple Resonance Probe	1070:1	 	

Resolution

Resolution is normally increased as the field strength of the magnet is increased. Overlapped peaks on the 300 NMR may be resolved at higher field strength. The partial ¹H spectra of a natural product are shown below. The spectra were collected on the 400 and 500 spectrometers within our NMR facility. The 800 Mz spectrum was collected at the University of Minnesota. The three peaks are completely resolved from each other at 800 MHz.



Dynamic Molecules

Molecules which are dynamic on the NMR time scale have a lower coalescence temperature on a lower field magnet. A single set of sharp lines is observed above the coalescence temperature. In the example below, the dynamic catenane compound had broad lines in the ¹³C spectra at room temperature. Due to sample constraints 50°C was the upper limit that the sample could be heated. On the 400 MHz NMR spectrometer, this temperature was not far enough above the coalescence temperature to result in sharp lines. However, considerably sharper lines were observed on the 300 MHz NMR spectrometer at 50°C due to the lower coalescence temperature.



IV. Software

Intro to VNMRJ

All of the computers in the facility are running VNMRj 2.2C and 2.2D. The interfaces of these two versions are nearly identical.

The program that communicates between the workstation and console is called the "acqproc". Occasionally this program loses communication between the console and workstation and needs to be re-started. Please ask facility personnel to do this if they are available. If no one is available, instructions are posted in the facility for re-setting the console.

Every group has a login account. They are responsible for remembering their user name and password. **Do not allow others to use the account or give out**

the password to anyone. If the facility determines that a user is allowing unauthorized users to use the account, privileges *will* be suspended.

The VNMR software currently uses linux for many things. Knowing a few simple linux commands can be very helpful (See Appendix for useful linux commands).

In a user account, there is a directory called "vnmrsys". This is a VNMR system directory. In this directory are many subdirectories. These directories are probably empty unless someone has put in a new macro, pulse sequence, shim file, parameter file, etc.

The most significant of the subdirectories:

psglib - pulse sequence library, contains uncompiled pulse sequences
seqlib - sequence library, contains the compiled, executable versions of pulse sequences
maclib - macro library
shims - stores personal saved shim files
parlib - parameter files
probes – probe files (created with an addprobe command)

shapelib – shaped pulse files (*filename*.RF)

gshimlib – gradient shimming maps and files (created from gmapsys)

Global directories are found in the **/vnmr** directory (a few directories up from the personal directories). The global directories, like **psglib**, contain all the pulse sequences, macros, parameters, shims etc. available to all users. To see the pulse sequences available on a specific machine, look in the **psglib** in the VNMR directory. The main VNMR subdirectories can be altered by facility staff only. Files may be copied into a personal **vnmrsys** and altered there. Any macros, pulse sequences, parameter files, etc. in a personal **vnmrsys** will be accessed preferentially to the global **vnmrsys** during operation.

To put a new pulse sequence into the **vnmrsys**, copy the uncompiled sequence (*seqencename*.c) into the **psglib**. Be sure to copy any necessary macros into the **maclib**, any parameters into the **parlib**, shaped pulses into **shapelib** etc.

To compile a sequence copied into the **psglib**,

>seqgen pulseq.c

This will compile the sequence and put the executable into the **seqlib**. Pulse sequences can also be compiled from the VNMR command line as **>seqgen(**'sequencename.c').

VNMR also has various packages for specific applications. The most common one in the facility is called BioPack and includes many pulse sequences for applications to biomolecules. These packages can be installed or activated in individual user accounts.

V. Data Files

Creating Data Directories

We ask that each group save its data within the data directory of the vnmrsys folder. The data directory is found at /home/usergroup/vnmrsys/data. The usergroup is designated by the group login. For example, if your group login is vnmr1, then the group data directory is found at /home/vnmr1/vnmrsys/data, but if your login in rhc, then the group data directory is found at

/home/rhc/vnmrsys/data. Please maintain individual directories for each member of the group. Directories can be made using the file browser. A file browser can be accessed under the Applications>System Tools.



After the file browser is opened, type the location of your group data directory in the location bar and create a new folder under the File menu.

Applications Places	System 🔟 🧒 🔆	r (* 19 10 10 10 10 10 10 10 10 10 10 10 10 10			🗐 12:08 PM 📢
	2	🕲 File Ejtit View Go	data - Fi <u>B</u> ookmarks <u>H</u> elp	le Browser	_ox
Computer	data on 129.15.12. 101	Back Forward	Up Stop	Reload Home	Computer Search
		Location: /ho	ome/vnmr1/vnmrsys/data		🐧 100% 🍳 🛛 View as Icons 🖨
vnmr1's Home		Places▼ ×			
		😴 Desktop	jon	matt	susan
	Data Station 1	S Floppy Drive			
VnmrJ		Gata on 129.15.12.1 Data on 129.15.13.3			
		Bata Station 1	05.9 GB		
Vomrl Admin					

The NMR facility automatically sends all files within the data directories of each group to a local server at the end of every day. Files that are not saved in the

data directory are not sent to the server. *However, the facility does not assume any responsibility for backing up user data files. Users must back up their own data.* Each spectrometer and data station computer has a flash drive. It is suggested that users transfer data via flash drives to their own computers where they can maintain a CD back-up of all data. Alternatively, files can be copied from spectrometer to workstation within the facility via ssh. See appendix A.6 for details in using the freeware Winscp to transfer data.

Creating a Desktop Shortcut to the Data Directory

It will be convenient to have a link on the desktop that will navigate directly to the data directory of your group.

Using the file browser, open the vnmrsys directory within the home directory and right click on the data folder. Click Make Link.



This creates a link inside the vnmrsys folder named link to data. Drag the link to data shortcut to the desktop.



Right click on the link to data and rename it to Data Directory.

ices	System 💼 🛞	®®\$\$
	2	Dpen Browse Folder
	data on 129.15 101	Open with "Cervisia" Open with Other <u>A</u> pplication
	Instructions_2	Cut Copy Capp Copy Capped Content Con
	Quata Station	Ma <u>k</u> e Link <u>R</u> ename
		Move to Trash
	Manual_08170 inprogress.do	Send to
	link to	Properties data

Creating the shortcut shown below:



Creating Data Transfer Shortcuts on the Red Hat Linux 5.1 Desktop



Click on Places

The Connect to Server box will open:

Set the Service type to SSH and fill in the Server, Folder and Name to user for connection

🎾 Co	nnect to Server X
Service type: Public FTP	
<u>S</u> erver:	
Optional information:	
<u>P</u> ort:	
<u>F</u> older:	
Name to use for connection:	
Browse <u>N</u> et	twork X Cancel Connect

Server-Internal IP	Server-External IP	Name for Connection
10.254.219.26	129.15.22.45	300 NMR
10.254.219.27	129.15.22.46	400 NMR
10.254.219.28	129.15.22.47	500 NMR
10.254.219.30	129.15.22.49	Data Station 1
10.254.219.31	129.15.22.50	Data Station 2

Note: Use the Internal IP address for the Server number if you are transferring between spectrometer and data stations. Use external IP address for the Server Number if you are transferring data to a computer outside of this network.

Folder: Location of your group's data directory

Generally: /home/login/vnmrsys/data

Specifically if your group login is rhc then the folder is /home/rhc/vnmrsys/data

Example: This is a shortcut made to Data Station 1 by the vnmr1 group.

	Con	nect to Server
Service <u>t</u> ype:	SSH	\$
<u>S</u> erver:		129.15.13.31
Optional in	formation:	
<u>P</u> ort:		
<u>F</u> older:		/home/vnmr1/vnmrsys/data
<u>U</u> ser Name	н	vnmr1
<u>N</u> ame to us	se for connection:	Data Station 1
🔀 <u>H</u> elp	Browse <u>N</u> etv	vork X Cancel Connect

The above information created the following shortcut:



Transferring the Data

To transfer the data double click on the desired short cut. This will open up box requiring your group password. Put in the password and click connect. Copy and paste your files from the local computer to the destination computer. *Example:* Transferring data from the 300 to Data Station 1 using the vnmr1 account: Click on the data directory shortcut on the desktop and navigate to the files that you wish to copy. NMR data have a .fid extension. These are actually folders containing four files. You must copy the entire .fid folder.



Open the appropriate folder to find the files that you wish to transfer.

System 🔞 🥱 🖏	S & S			📾 12:07 F
1		data - File Browser		;
<u>File E</u> dit <u>V</u> iew <u>G</u> o	<u>B</u> ookmarks <u>H</u> elp			
Back Forward	Up Stop	Reload Home	Computer Search	
Location: /h	ome/vnmr1/vnmrsys/data	\supset	Q 100% Q	View as Icons 🖨
Places▼ ×				
🔞 vnmr1				
🛷 Desktop	jon	matt	susan	

Click on the Data Station 1 shortcut to open a connection to the data station.



Enter in the password and click the connect button.

Applications Places	System 💼 🌏 🊳	\$880	2							12:08 PM 🚯
	٨	👏 <u>F</u> ile <u>E</u> dit	<u>V</u> iew <u>G</u> o J	<u>B</u> ookmarks <u>H</u>	data - Fil Ielp	e Browser				_ _ X
Computer	data on 129.15.12. 101	🖨 🗸 Back	Forward	ြာ Up	Stop	🤣 Reload	🙀 Home	D Computer	Search	1
		Loc	ation: /hom	ne/vnmr1/vnm	rsys/data			💐 100% 🍳	View as l	cons 🖨
vnmr1's Home		Places▼	×		7					
		😻 Desktop		jon		matt		susa	in	
		1	Authe	ntication Re	quired	×				
	Data Station 1		You must log	g in to access	vnmr1@12	9.15.13.31				
		<	Password:							
Vnmr] Admin			🗌 <u>R</u> ememb	oer password f	or this ses	sion				
			🗌 Save pa	ssword in <u>k</u> eyr	ing					
				×	Cancel	Co <u>n</u> nect				
CP4.lgz						1				

This will open a file browser on Data Station 1 (destination folder)

Applications Places System	• ™⊗⊗⊗ ≷{						🖮 12:31 PM 🔇
	1		dat	a - File Browser			_ _ X
	Eile	<u>E</u> dit <u>V</u> iew <u>G</u> o	<u>B</u> ookmarks <u>H</u> elp				
Computer data o	on 129.15.12.	ck <u>Forward</u>		Relead	Kome	D Computer	Search
		Location: /hor	ne/vnmr1/vnmrsys	/data		100% 🝳	View as Icons 🖨
vnmr1's Home		ces 🔹 🔪 vnmr1 Desktop File System	jon	ma	att	susa) n
8	data on 129.1	5.13.31 - File Bro	wser		. 🗆 🗙	l ocal	Folder
<u>File Edit View Go</u> Bookma	arks <u>H</u> elp					2004.	
Back Forward Up	Stop	🤣 🙀 Reload Hom	ی e Computer	Search			
Location: ssh://vnmr1	@129.15.13.31/ho	me/vnmr1/vnmrsys	/data 3 100%	🔍 🛛 View as Ico	ns 🖨		
Places▼ × iiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiii	jon	matt	susan				
Sile System	tion Folde	er					
lata on 129.15.1							
Data on 129.15.1 Data Station 1							
- 							
😻 🛛 🚳 data - File Browser		data on 129.15.13.3	31 - File Browser	💼 Starting Take S	Screenshot		9 🗖

Transfer the data from the local folder to the destination folder using the copy and paste functions under the edit menu. You may also "drag and drop" but be aware that you may be moving your files instead of copying them.

VI. Probe Files

Instrument Calibration

Proper operation of both chempack and biopack within the VNMRJ 2.2C and VNMRJ2.2D software requires the use of up-to-date and calibrated probe files. Setting up experiments from the software menu as described in the experimental instructions of this manual should result in the most recently calibrated values being automatically loaded into the experimental data sets. This can be checked by comparing the data parameter sets to the calibrated values within the probe files on the computer and/or the hard copies in the red calibration notebooks located next to the computer of each spectrometer.

Probe Files

Probe files are specific to each instrument probe. They contain all of the necessary calibration information. The actual calibration values can be found in the red calibration notebooks located next to each spectrometer and inside the text of the probe file itself. The location of the system probe files is in the /vnmr/probes/ directory of each instrument and can be viewed by any user. The system probe files can only be changed by the system administrator. Groups may manage their own local probe files in desired, but be aware that only the system probe files will be updated with current calibration data. Groups that create their own probe files will be responsible for their maintenance.

It is important for the user to be aware of which probe is on the instrument and which probe file is loaded in the software. The 400 and 500 spectrometers each have two probes. The news page of the faces on-line sign-up, and the first page of the red calibration notebooks next to each spectrometer have current probe installation information. The table below shows the name of the probe file that corresponds to each available probe in the laboratory.

Instrument	Probe	Probe File
Mercury VX-300	4-nuclei autoswitchable PFG probe	asw4405
VNMRS-400	Auto-X-Indirect Detection probe	AutoX_ID_8726
VNMRS-400	Auto-X-Dual Broadband probe	AutoX_DB_8790
VNMRS-500	Triple Resonance Probe	TR_8064 (for chempack use) HCN (for biopack use)

VNMRS-500	Indirect Detection Probe	id8297

Changing the Probe File

The loaded probe file can be seen and accessed at the bottom of the VNMRJ window.

Temp Spin Lock Probe 20.0 C 0.110 0.1 id8297	Idle Idle	Varian software VNMRJ VERSION 2.1 REVISION B.

Click on the Probe button to access the probe files.

			Probe		3
Sel	ect Prob	De		(System) ▼	1
	Tun	e sweep] 1	une gain 50]
	Edit Pr	obe			
	Edit	<u>U</u> ndo	Close	Abandon	

Click on the down arrow to view the probe file options. Select the desired file and close the box. 6

			Probe		2
Se	lect Prol)e			
	id8297			(System)	•
	HCN i d8297			(System) (System)	
	TR_8064 None	ļ		(System)	
	Edit Pr	obe			
	Edit	<u>U</u> ndo	Close	Abandor	ì
	Se	Select Prol id8297 HCN 1d8297 TR_8064 None Edit Pr	Select Probe	Probe Select Probe id8297 HCN id8297 TR_8064 None Edit Probe Edit Undo Close	Probe Select Probe id8297 (System) HCN (System) TR_8064 (System) None Edit Probe Edit Undo Close Abandor

VII. Gradient Shim Maps

All three spectrometers utilize gradient shimming. Gradient shim maps are made by the NMR facility staff and stored in the /vnmr/gshimlib/shimmaps/ directory. Users are welcome to make and store their own gradient shimmaps within their users gshimlib directory, but the details of making these maps is not found in this manual.

It is important that you are using a shimmap which corresponds to the probe installed on the instrument. You may also want to change shim maps to use the most current map. The name of the most current gradient shim map can be found on the instrument status report located in the three ring calibration notebook next to each of the spectrometer's host computers.

Loading a gradient shim map

The currently loaded gradient shimmap name can be found on the Standard Tab of the Acquire Page. The maps are named by the corresponding probe and the date that the map was made. In this case, the map name is

AutoX_DB_8790_lk_2011-02-25, meaning that it was generated on February 25, 2011 with the AutoX_DB_8790 probe on the magnet.

				VnmrJ		c
	<u>F</u> ile <u>E</u> dit <u>V</u> i	ew E <u>x</u> periments <u>A</u> cquisit	on Automation <u>P</u> roce	ess <u>T</u> ools <u>H</u> elp		
	: * 🖻 🗖	S O O			- L - M	
					From Concern	
	Set hardwar	e: operation complete			÷ 100	
	2C 🗾					
	Exp:3 S	eq: s2pul Index:	1			
\searrow	Er;					
\mathbf{i}	Ho					
	0					
		×			₽ ×	
×	Start Acqui	re Rocess Setup Hardware	ClearSampleInfo	Show Time Stop		
	Standard Lock	Operator : vnmi	1 Logout			
	Shim	Sample Name:		Sample Status	Find z0	
	Spin/Temp	SampleDir: (01)		Insert Eject	Gradient Autoshim	
		Lot Number:	Page:	Spin at 20 Hz	AutoX_DB_8790_Ik_2011-02-25	
	1	Notebook:	Eaddr.	Temp at 25.0 C	Autolock	
	-	Solvent DMSO		Spin: Off	lock Autoshim	
		Comment		Rate: 0 Hz	When: Never V	
	-	Automated pulse width		LOCK: NotReg	Shim method: z1z2 👻	
	<u>u</u>	Observe channel		Level:		

To change the gradient shimmap, first type gmapsys on the vnmrj command line.

. Coladadadadadadadadadadadadadadadadadadad		
<u>File Edit View Experiments A</u> cquisition	Automation Process Tools Help	
: 📩 🔛 🔜 🚱 🔘 🔘		- <u>1</u> - <u>1</u>
gmap sys		- 18aa
Gr gmap sys		
L Exp:5 Seq: s2pul Index: 1		

This will bring up the gradient shimmap pages. The option to change the map can be found on the Gradient Shim Tab of the Acquire Page.

<u> </u>	<u>le la </u>	<u> Alla Alla Alla Alla Alla Alla Alla All</u>			<u></u>			
e <u>E</u> dit <u>V</u> iew	E <u>x</u> periments <u>A</u> cquisiti	on Automation	<u>P</u> rocess <u>T</u> ools	<u>H</u> elp				
] 📂 🔝 🚳	00							Jr - W
mapsys								
Set hardware: of	peration complete							
*								
Exp:3 Seq: (gmapz Index:	1						
Start Acquire F	Show Time	Acquire St	op		Seq	uence]
Start Acquire F Defaults	Gridient Shim Setup	Acquire M	op lake Shimmap		Seq	uence Gradient Au	itoshim]
Start Acquire F Defaults Gradient Shim	Gridient Shim Setup	Acquire M	op lake Shimmap	bimman	Seq	uence Gradient Au Gradient	itoshim Autoshim a	
Start Acquire F Defaults Gradient Shim Acquisition Ulse Sequence Channels	Gridient Show Time Acquire Trial Spec	Acquire Store	op lake Shimmap Automake S	himmap	Seq	Gradient Au	itoshim Autoshim o	on Z
Start Acquire F Defaults Gradient Shim Poulse Sequence Channels Flags	Gridient Shim Setup Acquire Trial Spec Set Acquisition Parameter	Acquire Store	op lake Shimmap <u>Automake S</u> urrent mapname:	Shimmap Set by dat	Seq	Gradient Au Gradient Au	itoshim Autoshim o	on Z
Start Acquire P Defaults Gradient Shim Requisition Pulse Sequence Channels Flags Future Actions	Grdient Shim Setup Acquire Trial Spec Set Acquisition Parameter PFG H1 Homospoil	Acquire Store	op lake Shimmap Automake S urrent mapname: AutoX_DB_8790_Ik_	set by dat 2011-02-25	Seq	Gradient Au Gradient Au Gradient Quit Grad	itoshim Autoshim a dient Autosl	m Z
Start Acquire F Defaults Gradient Shim Acquisition Clase Sequence Channels Flags Future Actions	Grdient Shim Setup Acquire Trial Spect Set Acquisition Parameter PFG H1 Homospoil PFG H2 Homospoil	Acquire M tra s: iH1 iH2	Automake S Automake S urrent mapname: AuroX_DB_8790_Ik_ Load map AutoX_DB	bimmap Set by dat 2011-02-25 \$790_lk_2011	Seq	Gradient Au Gradient Au Quit Grad	itoshim Autoshim o dient Autosl	n Z
Start Acquire F Defaults Gradient Shim Acquisition Vulse Sequence Channels Flags Future Actions	Show Time Gridient Shim Setup Acquire Trial Spect Set Acquisition Parameter PFG H1 Homospoil PFG H2 Homospoil Gradient Type: nnc	Acquire Store	Automake S Automake S Automake S AutoX_DB.8790.Ik. Load map AutoX_DB Shims Used A	bimmap Set by dat 2011-02-25 .8790_Ik_2011 Window Size	Seq	Gradient Au Gradient Au Gradient Quit Grad Display Fit	toshim Autoshim o dient Autosl Display Shim Blot Shim	n Z
Start Acquire P Defaults Gradient Shim Togetiliter Vulse Sequence Channels Flags Tuture Actions	Show Time Gridient Shim Setup Acquire Trial Spec Set Acquisition Parameter PFG H1 Homospoil PFG H2 Homospoil Gradient Type: nnc Gabinn 21–24 first.	Acquire Store	Automako S Automako S Automako S Automako S Automako S Santa Santa Santa Automako S Automako S Auto	bimmap Set by dat 2011-02-25 8790_lk_2011 Window Size	Seq	Gradient Au Gradient Quit Grad Display Fit Plot Fit Set mappa	itoshim Autoshim o dient Autosl Display Shim Piot Shimn me into probe	n Z him nap
Start Acquire P Defaults Gradient Shim Topicius Pulse Sequence channels Tags Tuture Actions	Gridient Shim Setup Acquire Trial Spec Set Acquisition Parameter PFG H1 Homospoil Gradient Type: nnc Shim 21-24 first. Temp compensation of Temp	Acquire St M tra S: (H1 H2 (F F F F	Automake S Automake S Automake S Auto 208,8790.lk. Load map AutoX_DB Shims Used A	Set by dat Set by dat 2011-02-25 8790_lk_2011 Window Size Find Frequen	Seq	Gradient Au Gradient Au Quit Grad Display Fit Piot Fit Set mapna Add mapn	itoshim Autoshim o dient Autosl Display Shin Piot Shim me into probe ame into Gmag	n Z him map nap

Click the down arrow to see all of the available maps, highlight the desired map and it will load into the current mapname box.

Exp:3 Seq:	gmapz Index: 1		
tart Acquire F	Process Show Time Acquire	Stop	Sequence 🔍
efaults	Gradient Shim Setup	Make Shimmap	Gradient Autoshim
quisition Ilse Sequence	Acquire Trial Spectra	Automake Shimmap	Gradient Autoshim on Z
ags iture Actions	Set Acquisition Parameters: PFG H1 Homospoil H1	AutoX_DB_8790_Ik_2011-02-25	Quit Gradient Autoshim
	PFG H2 Homospoil H2 Gradient Type: nnc	Load map AutoX_D8_8790_Ik_2011 11-18-2010-lineshape # Shims Useq	Display Fit Display Shimmap
	Shim z1-z4 first.	ATgm ap Set WAutoX_DB_8790_Ik_2010-09	Set mapname into probe file
	Temp compensation off	Find Wind AutoX_D8_8790_Ik_2010-11 AutoX_D8_8790_Ik_2011-01	Add mapname into Gmap list
	Ignore spinner 🔹	Make Shimn AutoX_DB_8790_lk_2011-02 AutoX_DB_lk_2010-20-11	Reset Shims

Click the Quit Gradient AutoShim button to return to the regular VNMRJ pages.



VIII. Using the Varian Mercury VX-300 MHz Spectrometer

The 300 MHz NMR was purchased in 2000 under a NSF multi-user grant (CHE 0077707). It has a 4-nuclei autoswitchable PFG probe. It can collect 1H/19F/13C/P31 without retuning the probe or changing cables. This instrument is designed for walk-up use and is primarily used by the synthetic chemistry groups. Mornings, Evenings and Overnight time can be reserved in advance in large blocks, but daytime can only be reserved in 15 minute blocks.

300 NMR Sign-up Rules:

Mornings (8 am -12:00 pm)

Sign-up is allowed one day in advance for a maximum of 1 consecutive hour. Do not sign-up for one hour to use the instrument for 5 minutes. This time is reserved for longer experiments and inexperienced users.

Afternoons (12-noon to 6:00 pm) Sign-up is allowed only on the day of use and is restricted to a maximum of 15 consecutive minutes.

Evenings (6:00 pm to 8:00 pm) Sign-up is allowed one day in advance for a maximum of 2 consecutive hour

Overnight (8:00 pm to 8:00 am): Sign-up is anytime in advance for a maximum of 12 consecutive hours

Mercury 300 Tests and Assignments Certification

Student Name:

300-Test #1: Student gives demonstration of properly operating the instrument. The student must log into the spectrometer, insert sample, lock, shim, collect a 1H experiment, save the data and transfer the data to the data station. Speed is not required. At the completion of this test, the student is allowed to use the instrument Monday-Friday between the hours of 8 and Noon. Date Completed: Supervisor: Practical Assignment #1 Collect 2¹H NMR spectra of 2-Ethyl-1-indanone on the 300. This sample is located in the NMR laboratory. Please shim the sample as well as you can. Follow the instructions located in Appendix A.1 of this manual. Upon completion copy the assignment and give to the NMR facility staff. Date Completed: _____Supervisor: _____ Practical Assignment #2 Collect a ¹³C NMR spectrum of 2-Ethyl-1-indanone. Follow the instructions located in Appendix A.2 of this manual. Upon completion copy the assignment and give to the NMR facility staff. Date Completed: Supervisor: 300-Test #2: Student gives demonstration of properly operating the instrument. The student must log into the spectrometer, insert sample, lock, shim, collect a 1H experiment, save the data and transfer the data to the data station. The test must

Date Completed: _____Supervisor: _____

be completed within 15 minutes. Upon completion of this test, the student is

permitted to use the instrument at any time.

Collecting a 1D-proton on the Mercury VX-300 MHz NMR Spectrometer

Log in and double click on the VNMRJ icon.



Insert your sample clicking the Eject and Insert buttons on the Start> Standard Page.

andard ck im	Sample Notebook	Page	Insert Eject	Find Gradient Mapname:	z0 Shim	
in/Temp	Solvent	Acetone	Spin at 20 Hz 5	Lock Find P	asonanca	
	Std Proton pa	arameters		When:	Not used	Ŧ
			Lock	Shim or	Lock	
			Status: Regulated	When:	Not used	•
			Level: 27.7	Shim method:	z1z2	•

Set up a proton experiment by choosing proton under the experiment menu. If there are no experiments under the Experiment menu, uncheck the Use Study Queue option under the Acquisition menu.

La	Vnm	nrj
<u>F</u> ile <u>E</u> dit <u>V</u> iev	Experiments Acquisition Automation Process	<u>T</u> ools <u>H</u> elp
: *) 📂 🗔 🧕	Setup NEW parameters for	
	Proton	
	Carbon	
ArrayedSpectra	Fluorine	
1D 2D	Phosphorus	
Viewport	Other Nucleus	•
Frame	H1 Relax, Measurements	•
[Holding		
₽ ×	Convert current parameters to do	
Std1D (HC)He	Standard 1D experiments	•
Jn(CH)corr	Solvent Suppression – Select peaks	•
(HH)Homo2D	📕 Homonuclear Correlation Experiments	•
Common	Indirect Het. Corr. (Basic)	•
	Indirect Het. Corr. (More)	•

Start Acqui	re Process Setup Hardware Show Time	Stop	
Standard Lock	Operator mta	Insert Eject	Find z0 Gradient S
Shim Spin/Temp	Notebook Page	Spin at 20 Hz 🖌	Mapname:)B_87
	Solvent Chloroform DMSO D2O other D2O	Temp at 26.0 C	Lock Find Res
	Comment	Lock	Shim on L
	STANDARD PROTON PARAMETERS	Status: NotReg	When: N
		Level:	Shim method: z

Select the desired solvent on the standard tab of the Start page.

Lock:

Click the Find z0 button on the Lock tab of the Start page.

	A	В	C	
Start Acquire Process	Setup Hardware Show Time	Stop		Ф ×
Standard Lock On	Off Insert Eject	Lock Scan	Spin On Spin: Off Spin Off 0 Hz	
Shim Status: Not	Reg ZO	_		
Spin, remp	-12314 ±1 +12314			
Lock Level:	32 ±1 32			_
	27 ±1 27			_
	Phase ±1 269			_
Find z0	Gradient Shim	Set Up Hardware	Load Shims Into Hardwa	re

Click the Lock Scan, Lock Off and the Spin on buttons. Adjust (raise or lower as necessary) the Power and Gain until the signal completely appears on the screen.



Adjust z0 to maximize the lock signal. For best results adjust z0 until signal is maximized, then adjust phase until the signal is maximized. Clicking the middle

wheel of the mouse will change the steps of z0, power, gain and phase to 1, 10 or 100. If the lock level is fluctuating at a low power level, then lower the gain.

Eyp:2 Seq: Proton Index: 1
Start Acquire Process Serup Hardware Show Time Step
Skindard Lock the Off Insert Eject Lock Scan Spin On Spin. NotReg
Shim Status Regulated Spin/Temp 20 +1 392
Jack Levels
Phase at the second sec
195 11 135
Find zo Gradient Shim Load Shims Into Hardware
×
teractive Set hardware: operation complete

Turn the Lock On when you are finished locking. <u>You must</u> unclick the Lock Scan button for the other buttons to become available.

Shim by clicking the Gradient Shim button on the Lock tab of the Start page. If you want to manually shim z1-z5, click the shim tab. It should not be necessary to shim the higher order shims.

Start Acqu	uire Process Setup Ha	rdware Show Time Stop	Ф ×
Standard Lock Shim Spin/Temp	Lock On Off Status: Regulated Lock Level: 46.4	Insert Eject Lock Scan Spin On Spin Off O Hz 20 -9397 -	
	Find 20	Gradient Shim Set Up Hardware Load Shims into Hardware	
90	ldie	Gradient autoshimming on Z done! 1 iteration	

Alternatively, you may manually shim using the buttons on the Start>Shim Page to maximize the lock signal.

andard Lock Scan		21	± 10	×1	± 1 XY	±
im im in (Tame 20		363 Z2	± 1	<u>-173</u> Y1	351 ±1	+
392 ±1	40.4	326 Z1C	+ 1 23		<u>262</u>	
Lk Power 12 ±1		319	154	± 10 99		2
Lk Gain		300	±1	± 10	± 1 _45	±
<u>19</u> ± 1			Z5 -114	± 10		
Lk Phase ±1	Lock 4	10.4 Sni	n On Spin Regulated	Recr gain		

Start your experiment by clicking the Acquire button. Parameters may be accessed under the Acquire pages. Basic parameters such as sweep width, number of scans, relaxation delay, etc are found in the default 1H tab. Automatic plotting and integration can also be set to on or off here.

A		B	c
Start Acquire Process Acqu	re Acquire & Transform	Stop Show Time Sequence	Arrays
Default H1 Acquisition Pulse Sequence Channels Flags Future Actions Relaxation Delay Number of Scans Spin Tuge method	om) $14 \rightarrow -2$ ▼ 0 Upfield -2.0 ▼ ees] 45 ▼ anje 45 ▼ 8 ▼ 20 Hz lohi ▼	Transform size: 32k Line Broadening [Hz] Plot when done Spectrum: Full Parameters: Basic, Top Le Peak Values: Integrais: Partial	

An "old style" parameter page can be accessed by typing dg in the command line and clicking on the Process>Text Output Page



More instructions regarding changing the parameters of a basic 1D experiment can be found in section X of the this manual: Changing Basic 1D Parameters.

Save your data by choosing File > Save As from the drop down menu.



This opens a box in which you can type your filename. Clicking the Home icon will take you to the data directory within your account. You can create folders here using the create folder icon. Double click on your folder.

	- Lunterd -		Save				
Choose Home Dir	ectory/home/m	nta/vnmrsvs/dat	a				
Dir 1		Dir 2		Dir	3	Dir 4	
Save Inc. 🗖 dat					- (-
Save In.	d						
📑 susan							
inpstudy.							
						/	
					/		
File Name:				_ /			
Files of Type:	fid						
These Types				/			
						Save	Cance
				12	1	ó i	
			/				
			/				
		•					
one				1	Crea	te New	
		1 A	l cat	4	Folde	<u>ə</u> r	
ry		╵┃╘┳┨			1 0100	51	
				٦			
	G	ο Το Ηοι	me				
				1			
	Di	rectorv					

Eject your sample clicking the Eject and Insert buttons on the Start> Standard Page.

andard ock	Sample			Insert	Fjert	Find	z0	
nim pin/Temp	Notebook Solvent	Acetone	Page	Spin	at 20 Hz 🗹	Mapname:	Shim	
	Comment			Temp	at 3.0 C 🗌	Lock Find R	esonance	
	Std Proton pa	rameters				When:	Not used	•
				Lock		Shim or	n Lock	
				Status:	Regulated	When:	Not used	•
				Level:	27.7	Shim method:	z1z2	•

Close VNMRJ software

Type exit in the command line to close the VnmrJ software.



Log off Computer

Click Log Out from the Systems Menu.



IX. Using the Varian VNMRS-400 MHz Spectrometer

The Varian VNMRS-400 NMR Spectrometer was purchased in 2007 by the multiuser NSF grant CHE#0639199. It is equipped with two probes: Auto-X indirect detection probe and Auto-X dual broad band probe. Please request a probe change at least 2 working days in advance.

Training will begin on the 400 upon successful completion of 300 training. Certification to use the 400 NMR will be given as NMR Test #400-1, and Practical Assignments #400-1 and #400-2 are completed as approved by Dr. Nimmo.

400 NMR Sign-up Rules:

Rules for Advanced Signup

Advanced signup is defined as signing up for time more than 24 hours in advance. These rules do not apply if this time has not been taken within 24 hours of the start of the experiment.

Weekdays (Monday – Friday 8:00-4:00):

Advanced sign-up of this time is reserved for experiments which must be monitored, 400 beginning students who are only allowed to use the instrument during this time and experiments which require help from the NMR facility staff.

Weekends: (Friday 6 pm – Monday 8 AM) Please limit advanced sign-up to 24 hours during the weekend unless approved by the NMR facility staff.

400 NMR Tests and Assignment Certification

Student Name:

400-Test #1:

Student must give demonstration of properly operating the instrument including sending the data to the data station. This test includes operation of the robot, automated tuning unit and shimming the standard line shape sample to instrument specifications.

Date Completed: _____Supervisor: _____

400- Practical Assignment # 1 Measure the ¹H pw90 for any sample. Follow the instructions located in Appendix A.3 of this manual. Make a copy of the assignment and hand it in to the NMR facility staff.

Date Completed: _____Supervisor: _____

400- Practical Assignment # 2

Collect a gCOSY experiment using any sample. Follow the instructions located in Appendix A.4 of this manual. Make a copy of the assignment and hand it in to the NMR facility staff.

Date Completed: _____Supervisor: _____

Upon completion of 400- test #1, practical assignment #1 and practical assignment #2, the student is permitted to use the instrument at any time.

Collecting a 1D-proton on the Varian VNMRS-400 MHz NMR Spectrometer

Put your sample tube in a spinner turbine and set the depth using the depth gauge. Put the sample into one of the numbered slots in the sample tray and take note of the slot number. **Do not use slot #1 (upper left corner) or slot #0 (lower right corner).** A standard sample should be in the magnet which will be taken out and put into slot #1 when your sample is inserted. Slot #0 must be left empty. If the robot takes your sample out and places it in slot #0, you must remove it.



In this example, the sample is placed in slot #19.

Login and start VNMRJ by left double clicking on the VnmrJ icon on the desktop.



Set up a proton experiment by choosing proton under the experiment menu. If there are no experiments under the Experiment menu, uncheck the Use Study Queue option under the Acquisition menu.



Select the desired solvent on the standard tab of the Start page.

Standard Lock	Operator mta	Insert Eject	Find z0 Gradient Sh
Shim Spin/Temp	Sample Notebook Page Solvent Chloroform DMSO D20 other D20	Spin at 20 Hz V Temp at 26.0 C	Mapname:)B_879 Lock Find Reso When: Sin Shim on Lo
		Status: NotReg Level:	When: No Shim method: 21

Insert your sample by typing loc=slot# change. In this example the sample was placed into slot #19, so loc=19 change will cause the robot to eject the standard sample, replace it into slot #1 and then insert the sample from slot #19 into the magnet.


Lock:

Click the Find z0 button on the Lock tab of the Start page.

	A		В	C
Start Acqu	ire Process Setup H	ardware Show Time 🧾	Stop	Ф X
Standard Lock	Lock On Off	Insert Eject	Lock Scan	Spin On Spin: Off
Shim Spin/Temp	Status: NotReg	²⁰ ±1 .12314 ==		
	Lock Level:	$\begin{array}{c c} \hline & -12514 \\ \hline & Power \\ \hline & 32 \\ \hline & 52 \\ $	(
		27 ±1 27 = Phase ±1 269 =		
	Find z0	Gradient Shim	Set Up Hardware	Load Shims Into Hardware

Click the Lock Scan, Lock Off and the Spin on buttons. Lower the Power and Gain until the signal completely appears on the screen.



Adjust z0 and Phase to maximize the lock signal. The lock signal is gold and the phase difference is blue. For best results adjust z0 until signal is maximized, then adjust phase until the signal is maximized. Continue alternating between z0 and phase until the signal doesn't change when the Lock is turned on. The phase difference (blue) should be zero (flat line) when the lock signal (gold) is maximized. Clicking the middle wheel of the mouse will change the steps of z0, power, gain and phase to 1, 10 or 100.

et hardwa	are: operation comp	ete	
Exp:1	Seq: Proton	Index: 1	
	/		
tart Acqu	uire Process Setur	Hardware Show Time Stop	ņ
andard ock	Lock On Off	Insert Eject Lock Scan Spin On Spin Off	Spin: Regulated 20 Hz
oin/Temp	Status: Regulati	a Z0	
	Lock Level:	rower ±1 39	
	57.0	33 +1 22 Phase ±1 305	

Turn the Lock On when you are finished locking. You must unclick the Lock Scan button for the other buttons to become available.

Shim by clicking the Gradient Shim button on the Lock tab of the Start page. If you want to manually shim z1-z5, click the shim tab. It should not be necessary to shim the higher order shims.

Start Acquire	Process Setup Ha	rdware Show Time Stop	з Х
Standard Lock	ck On Off	Insert Eject Lock Scan Spin On Spin: Off O Hz	
Shim Spin/Temp L	tatus: Regulated .ock Level: 46.4	Z0 -9397 9397 - Power - 3a - Gain ±1 33 + Phase -	
	Find z0	Gradient Shim Set Up Hardware Load Shims Into Hardware	x
790	idie	Gradient autoshimming on Z done! 1 iteration	

Tune your sample by typing protune in the command line. (Alternately you may use the protune('calibrate') macro: see page 42.)

v	VnmrJ
<u>File Edit View Experiments A</u> cquisition Pro	ocess <u>T</u> ools <u>H</u> elp
: 📩 🖻 🔜 🚱 🔘 🔘	
2D Cryo ArrayedSpectra	Gradient autoshimming on Z done! 1 iteration
Protocols Frame Viewport 10	- protune
Experiment Panel	
Std 1D Home 2D History 2D / Sel 1D	o ⊏xp.i Seq: Proton Index: 1
Std ID Hollio 2D Hetero 2D Ser ID	

Start by tuning the high band to H1 (or F19). If tuning proton, then type H1 in the nucleus box under Advanced Tune and choose Fine for the Tune Criterion. Click Tune to Criterion to begin the tuning.

📤 Applications Actions 🥫 🔗 🥸 🍣 🌄 資	
 Tune Probe 	×
Quick Tune P31 C13 N15 P31 C13 N15 P31 C13 N15 P31 C13 N15 P31 C13 N15	
Tune to Criterion	
\sim \sim	
Edit Undo Close Abandon	

The acquisition status will return to idle and a Tuning done ok message will appear when tuning is completed.

]	Find z0	Gradient Shim	Set Up Hardware	Load Shims Into Hardware	×
79	ldle	Tuning done:	ok – tuned to 399.967 Mhz w	ith match at within 0.1 percent of optimum	1

You may also tune the low band at this time (i.e. C13, P31, etc) by the same method described above. Click the close button at the bottom of the window when finished.



Start your experiment by clicking the Acquire button or typing ga in the command line. Parameters may be accessed under the Acquire pages. Basic parameters such as sweep width, number of scans, relaxation delay, etc are found in the default 1H tab. Automatic plotting and integration can also be set to on or off here.

Image: Spectral Width [ppm] 142 Transformation Acquisition Puise Sequence Channels Flags Future Actions Spectral Width [ppm] 142 Unewfordening [Hz] Prise Sequence Channels Future Actions Downfield [14.0) Upfield P2.0 Unewfordening [Hz] V Puise Width (degrees) 45 V Plot When done Spectrum: Future Actions Relaxation Delay (Sec) 1 V Parameters: Easic, Top Left V Number of Scans 8 V Peak Values: V Tune method Ibhi V Partial V	Start Acquire	Process Acquire	Acquire & Transform	Stop Show Time	Sequence Arrays	T X
	Default H1 Acquisition Pulse Sequence Channels Flags Future Actions	Proton Spectral Width (ppm) Downfield [14.0 Puse Width (degrees) Enter pulse angle Relaxation Delay (sec) Number of Scans Spin Tune method	14 → -2 ↓ Upried [-2.0 45 ↓ 45 ↓ 8 ↓ 20 Hz John ↓	Transform vice: Lins Hoadening (H Plot when done Spectrum: Parameters: Peak Values: Integrals:	z] Full v Basic, Top Left v Partial v	

Save your data by choosing File > Save As from the drop down menu.



This opens a box in which you can type your filename.

✓	Si	ave		3
hoose Home Directory: /home	/mta/vnmrsys/data			
Dir 1	Dir 2	Dir 3	Dir 4	
Save In: 🗖 data			- a d d <mark>8</mark>	3=
🗂 susan				٦
🗂 tmpstudy				
File Name:				
Filer of Type: fid				5
Thes of Type. Ind				<u> </u>
			Save Cancel	
				_

To eject your sample, type loc=1 change in the command line. This will eject your sample, replace it into its slot (slot#19 in this example) and put the standard sample from slot #1 into the magnet.

-		2						
2	2D Стуо	Ar	rayedSpectra		/	ds		
	Protocols		Frame	Viewport	1D	▼ loc = 1 d	change	
Đ	xperiment Pa	nel		K.	Φ×	o Exp:1	Seg: Proton	Index: 1
Ĺ	Std 1D Ho	mo	2D Hetero 2D	Sel 1D		C LAPIT		

Type exit in the command line to close the VnmrJ software.

		2					
2D	Cryo	ArrayedSpectra		/	ds		
Ĺ	Protocols	Frame	Viewport	1D	exit)	
Expe	riment Pa	inel		X		Soguinfo	Index: 1
Std	1D Ho	mo 2D Hetero 20	Sel 1D		U LAP.I	Seq. IIII0	IIIGEX. 1

To completely log off, choose Log Out from the Actions Menu.



Tuning the Broad-Band Probe Using protune('calibrate')

The protune macro will automatically tune 1H/13C/19F/31P. However, the broadband probe is tunable to many additional nuclei. In this example, we will be tuning the probe to ⁷Li which is at 155.440 MHz on the 400 MHz NMR.

Type protune('calibrate') in the command window of VNMRJ.



This will open the manual calibrating program as shown below.



Type the desired frequency in the 'Tune to' box and then click the 'Tune to' button. (155.440 for this example).



The tuning will begin and continue until the bottom of the frequency dip reaches the bulls-eye.



Click the Quit button when the tuning is finished.

X. Using the Varian VNMRS-500 MHz Spectrometer

The 500 NMR spectrometer was originally purchased as a VXR model in 1987 and was completed upgraded to a VNMRS model by the University in 2006 with the addition of Dr. Susan Schroeder and Dr. Robert Cichewicz to the Department of Chemistry and Biochemistry Faculty.

The spectrometer is a three channel, 28 shims Varian VNMRS-500 equipped with two probes: triplet resonance $H^{13}C/^{15}N$ PFG probe, and a tunable indirect detection ${}^{1}H^{15}N-{}^{31}P$ PFG probe.

Training will begin on the 500 upon successful completion of 400 training and Test #500-1. Certification to use the 500 NMR between 8-5 Monday-Friday will be given as NMR Test #500-2 is completed as approved by Dr. Nimmo. Unlimited access will be given upon the completion of NMR Test #500-3

Advanced sign-up rules

Advanced sign-up for less than two hours will be limited to Wednesday afternoons. No one person is allowed to sign up *in advance* for more than 18 consecutive hours, with this exception of Dr. Schroeder's group which has ownership of 50% of all the available NMR time.

Operation Instructions

The 500 NMR is operated using VNMRJ2.2C software which is identical to the 400 NMR. With the exception of the autosampler and tuning, operation of the 500 NMR is consistent with the instructions presented in the 400 NMR operation section.

VNMRS 500 Tests and Certification

Student Name: _____

500-Test #1:

The student will be given a written test given by Dr. Nimmo. This test will be similar to the test at the end of the appendix. Upon the successful completion of this written test the practical training sessions will begin.

Date Completed: _____Supervisor: _____

500-Test #2

The student will demonstrate the proper use of the instrument. This will include setting the depth of the sample and inserting the sample into the magnet, changing the probe file, setting up a 1H experiment and tuning the probe. The student will be expected to properly identify all probe cables and connections and equipment components. At the successful completion of this test the student will be allowed to use the NMR from 8am – 5 pm, Monday –Friday.

Date Completed: _____Supervisor: _____

500-Test #3

The student will demonstrate the proper use of the instrument to a NMR faculty member. This will include setting the depth of the sample and inserting the sample into the magnet, changing the probe file, setting up an 1H experiment and tuning the probe. The student will be expected to properly identify all probe cables and connections and equipment components. This will also require the signature of the student's research director. At the successful completion of this test the student will be allowed to use the NMR at any time.

Date :	NMR Faculty:
Date:	Research Director:

Tuning the VNMRS-500 MHz NMR probes

WARNING: If you are a new user and you don't know how to do this or are unsure, **GET HELP**. Do not attempt to just "figure this out". Equipment can be damaged. Tuning probes may not be allowed in other NMR labs so please check with the lab manager before attempting to tune probes elsewhere.

The VNMRS model has some differences to earlier models. These differences result in some differences in tuning of this instrument and other models of instruments. There is no need to move any cables for the tuning process. Tuning the probe reduces the reflected power. Probe tuning requires adjusting the tune wands very carefully to reduce the reflected power to the 50 range on the tune interface box. Tuning the probe will insure that the 90-degree pulse width determination will be accurate and the pulse width will be as short as possible. If the sample is high in salt, the probe may not tune down to 50.

Tuning Tips: The wands can be very sensitive. Try to make small adjustments. If minimizing both the tune and the match and the reading doesn't go any lower, try moving the match part either clockwise or counter-clockwise to slightly increase the signal. Then try to minimize the signal again using the tune part. If the signal improves, continue to move the tune part in the same direction and repeat. If the signal does not improve, try adjusting the tune part slightly in the opposite direction to increase the signal and try to minimize again with the match part. This process should helps get out of any local minimum it is trapped in. Do not force the tuning rods, if they feel "stuck" stop and ask for help.

To tune the probe, first setup the correct frequencies through the computer. To do this, load the experimental parameters, i.e. setup a proton experiment then type **su** on the command line. Look at the channels page of the Acquire tab to see which nuclei is set up on each channel. Channel 1 is 1H, and channel 2 is carbon.

Start Acquire F	Process Show	Time Acquire	Stop MoveSW	Est pw90	Sequence RefreshPanel
Default H1 Acquisition Pulse Sequence Channels Flags Future Actions Overview	PROTON Channels: Nucleus / Freq. Offset Dec On/Off Dec Modulation 90 Degree at Pwr = dmf Waveform at resolution	Observe 1 H1 +99.882 MHz H99.9 Hz Homo per 1 offset 0.0 Hz × S.70 us at	Display Sequence Decouple 2 C13 125.707 HHZ 0.0 HZ 1500.0 us at 0 200 parp.1 1.00 degrees	Arrays Decouple 3 MHz 0.0 Hz • n c 31.0 us at 1 = \$2258 1.00 degrees	
	Ref. sample	Me4Si, 1% in CDCI3	Me4Si, 1% in CDCI3		

Set the tuning box to channel 1 and check that the attenuation is set at 9. If the meter reading is very large or off-scale, lower the attenuation. Tune the proton wand (top part is tune, bottom is match).



Adjust the wands to minimize the value of the reflectance power. Raise the attenuation back to nine (if necessary) as the reading decreases. Users should be able to reach a reading of 50. Always end with the attenuation at 9. Change the channel to 2 and tune the carbon wand if desired. When you are finished the tuning box must read channel 0 attenuation 9. The tuning box is only accurate to a value of 50; it may be helpful to use mtune to achieve a better tune.

Using mtune

To start the mtune program, type **mtune** on the command line. To observe the high band tune choose Tune RF channel 1; center frequency H1 and click Start Probe Tune.



To change to low band, click Stop Probe Tune, change Tune RF Chanel to 2 and center frequency to desired nucleus (C13 in this case). Click Start Probe Tune. Chanel 3 can be tuned to 15N in the same manner.



XI. Changing Basic 1D Parameters

A list of parameters can be accessed on the text output page of the process panel after typing dg in the command window. A short description of each parameter listed can by found in the Command and Parameter Reference Manual. Each parameter can be changed by typing the abbreviation = value in the command window. For example to change the number of scans from 1 to 16, simply type nt=16 in the command window. It is necessary to type dg for the list to refresh, showing that the change has been made.

Selected	a parameter list
at	Acquisition time
d1	Delay time between scans in seconds –recycle delay
ni	Number of increments
np	Number of points
nt	Number of transients
pw	Pulse width
pw90	90° pulse width
seqfil	Pulse sequence being used

Selected parameter list

sfrq	Observed frequency (frequency of tn)
SW	Sweep width (of tn)
tn	Transmitting nuclei (observed nuclei)
tof	Tuner offset frequency—center of the spectrum
tpwr	Transmitting nuclei power level
dfrq	Decoupling frequency
dmf	Decoupling Modulation Frequency
dn	Decoupling Nuclei
dof	Decoupling Offset Frequency
dpwr	Decoupler Power

Setting the frequency window (adjusting sw and tof)

The transmitter offset (**tof**) will always be in the *center* of the spectral window (**sw**). The spectral width (**sw**) of a spectrum may need to be altered and with it, the transmitter offset (**tof**). Use the smallest spectral width that accommodates the expected chemical shifts of the sample. Beware of using one too small though because peaks that are cut off will result in folded peaks.

When using water samples, set **tof** on the water peak. To do this, set the cursor on the water peak and type **movetof**. To move the spectral window, enclose the spectrum with the cursors, expand, and type **movesw**. This command moves the **tof** to the *center* of the new spectral window so **tof** will change. Manually reset the **sw** by using **>sw=#** if the **tof** needs to stay the same. The **sw** should be held constant for a series of experiments.

To select a spectral window, decide the approximate range where the signals will fall and select a slightly larger window. If all the signals can be observed in one scan, the window can be set manually using the **movesw** command. The **sw** is in Hertz. ppm = sw/sfrq. Therefore on the 300 NMR spectrometer, **sw=3000** is a window of approximately 10 ppm for a 1H spectrum. For a carbon spectrum on the 300 NMR spectrometer (carbon **sfrq** ~ 75 MHz), **sw=3000**, **sfrq=75**, is a window of approximately 40 ppm. Remember the spectral window is not a specific frequency range, so **sw=3000** could be 0-10ppm or -1 to 9 ppm or 1-11 ppm etc. The specific frequency range is set by the location of the transmitter offset (tof). Tables of tof values are found in the red calibration notebooks next to each spectrometer.

Changing the sweepwidth using the **movesw** command:

Simply place the cursors around the peaks and type movesw. This will automatically change the sw to the position of the cursors and adjust tof to the half-way point between the cursors. In the spectrum below, the sweepwidth was changed to 0 - 9 ppm with the center at 4.5 ppm. This corresponds to a sw of 2700 Hz and a tof of -184.9 on the Mercury VX-300 NMR Spectrometer.



Changing the sweepwidth and tof manually:

The transmitter offset (**tof**) will always be in the center of the spectral window (**sw**). The tof is a distance in Hertz, plus or minus, from the original center and will not match the scale when swapped from ppm to Hz. The tof parameter is spectrometer dependent and will therefore vary from instrument to instrument. Therefore, two different 500 NMR spectrometers will have different tof values. Tables of values for tof values with their corresponding ppm values are found in the red calibration notebooks next to each spectrometer.

To manually set tof, place the cursor at the desired location of the spectrum center and type movetof. Then type sw=#(hertz) for the desired window size. For example, if I wanted to have a 10 ppm window on the 300, with the center at 8 ppm, then I would set the cursor at 8 ppm and type movetof sw=3000. I would then collect the spectrum and reference the peaks. Alternatively, I could look at the tof value in the calibration notebook and type tof=# sw=# on the command line.

Setting the first delay (d1) and acquisition time (at).

The delay time, **d1**, is the time between the end of the acquisition time, **at**, and the next element in the pulse sequence. The **total recycle time** is the time between the last pulse and the next pulse at the beginning of the pulse sequence, typically is **d1 + at**. The total recycle time usually reflects the T1, longitudinal relaxation time, of the sample. Typically, a value of 1-2 seconds will suffice. If the T1 is unknown, a quick T1 experiment should be done to estimate it. The T1 can be long especially in cases of heteronucleus detection such as carbon or other X nucleus. T1 values for quaternary carbons can be very long and the total recycle time should reflect this.

The acquisition time, **at**, is usually set based on a rough approximation of T2, transverse relaxation time. The T2 indicates how long the signal in the FID will last. As a general rule, large molecules such as proteins relax very quickly and require a short **at** (such as **at=0.5**). Small molecules relax much more slowly and the FID can last several seconds so often using a longer **at** is favorable (such as **at=3 to 10** seconds). If you type dps on the command line the pulse sequence will be displayed. This is a graphical representation of the pulses that are taking place during each scan.



To change d1 or at, type d1=# and at=# on the command line. Type dg to refresh the parameter list and dps to refresh the pulse sequence.

XII. Basic 1D Processing

Opening Saved Data

Click Open..(under File in the Main Menu).



This will open a pop-up window. Clicking Home will take you to the data directory within your account. You can create folders here using the create folder icon. Double click on your folder.

w Experiments Acquisition Au	itomation Process Tool	ls Help		
	0	pen 🥂	×	
Choose Home Directory: /home/v	/nmr1/vnmrsys/data		•	Go To Home
Dir 1	Dir 2	Dir 3	Dir 4	Directory
Look In: 🗖 data			▼ 🖬 🗂 🗂 🔡 💳	
calibration_data			Home	
Danny_moore				
murali			Golupiono	- JA 🔂 🗖
Tanui 🗖 rice			directory	
🗖 susan			directory	7
				/
				Create New
				Folder
File Name: HP-20-MEOH-2	-HSQCAS-400MHZ.fid			
Files of Type: All Files			-	Highlight your
			Open Cancel	spectrum and click Open.

Spectra have a .fid extension.

			Open	in in the second se
Choose Home D)irectory:/home/vnm	nr1/vnmrsys/data		-
Dir	1	Dir 2	Dir 3	Dir 4
Look <u>I</u> n: s	usan			- A C C 8 =
Shekar_40	0_h1.fid 00.fid 00_t2.fid 20407_bc.fid fid			
File <u>N</u> ame:	styrnine_400.fid			/
Files of <u>T</u> ype:	All Files			
ŝt				Open Cancel

The spectrum will processes automatically and appear on the screen. Use the toolbar icons to access cursors. If the spectrum seems to "freeze", either click the redraw icon on the toolbar or type ds in the command window. The vertical scale can be changed with the mouse wheel or by typing in a numerical value for vertical scale (i.e. vs=5000).



Toolbar Icons

•

Ŧ

n

λL	Access one cursor, click again to access a second cursor
	Show the full spectrum
X	Reset to full display
	Zoom In
<u>_</u>	Zoom Out
٩,	Zoom More
<u></u>	Pan and Stretch More
<u>,,</u>	Integration
<u></u>	Show/Hide Scale
JUL	Show/Hide Threshold
16	Phase Mode
2	Redraw Spectrum
9	

Phasing

To display the spectrum on the screen, type **ds** Type **f full** to see the full spectrum in the full screen. The spectrum will likely need to be phased. To phase the spectrum, click on the

phase icon on the toolbar **M**. Use the left mouse button to click on a portion of the spectrum. Hold down the left mouse button and move the mouse to flatten the spectrum's baseline and makes the peaks symmetrical. The right mouse button will make fine adjustments and the center mouse button will increase or decrease the vertical scale. Simple spectra can be phased crudely with by typing **aph** (automatic phasing).

The spectrum below is a newly acquired spectrum which needs to be phased.



This spectrum has been "wrapped around" by bad phasing. It can be corrected by setting lp=0 and rp=0 and then phasing again.



Beware of becoming dependent on **aph** since autophasing does not always work and often leaves small errors in the phase that are better corrected manually. The **Ip**, or left phase, value corresponds to a first order phase correction that is frequency dependent, meaning it affects peaks differently. The **rp**, or right phase, value corresponds to the zero-order phase correction and is an error seen in all peaks.

Referencing the Spectrum

The spectrum may be manually referenced using the cursor. Expand the region around the peak you wish to reference by placing a cursor on either side of the peaks and clicking the zoom in icon.



Lower the vertical scale to bring the desired peak to scale. This can be done by clicking the mouse wheel under the scale (To raise the vertical scale click the mouse wheel at the top of the black part of the screen).





Put the cursor on the peak. In this example, it is the chloroform peak.

Type in the reference number on the Display page of the Process Tab.



Weighting Functions

To add a weighting function to a FID, type **wti** on the command line. This opens the weighting function screen. The mouse controls are shown on the screen. The right mouse button toggles the spectrum (top box) off and on, the center button controls the vertical scale of the FID (bottom box) and the spectrum, and the left mouse button controls the adjustment of the weighting function (green line in the center box).

Start by clicking on line broadening icon on the toolbar. This should bring up an exponential function in the middle window (green line). Click on the right mouse button in the top box to see the spectrum. Use the left mouse button to adjust the green line (weighting function) and observe the change in the spectrum. The value of the line broadening (**Ib**) is shown on the bottom of the screen.

pe Referei	nce frequency 7.2400	02 ppm, (2895.73 Hz)					
Exp:1	Seq: PROTON	Index: 1					
						Mul	
	~~~~~						3
Mouse bu Vf	ttons: Left - weighting, Cer vs 2:1b 3057.979 0.675	nter – vf/vs, Right – spectrum Bisb 4:sbs	on/off 5:gf 6:gfs	7:awc			
	0.075	unusea unusea	unusea unusea	0.004	_		<b>#</b> ×

The interactive weighting screen displays the FID, spectrum, and the weighting function applied. Above, an exponential function, line broadening (**Ib**), is applied to the FID. Line broadening improves signal-to-noise at the expense of resolution. Larger values also can improve FID truncation artifacts. Negative values of Ib can improve resolution, compared to the unweighted spectrum, at the expense of signal-to-noise.

After a good adjustment is found, type **wft** in the command line. This performs a weighted Fourier transform and applies the weighting function. The spectrum may need to re-phased slightly after this.

A sinebell function can also be used for resolution enhancement. To use a weighted function like a shifted sinebell, go to the interactive weighting screen by typing wti on the command line. Turn off any weighting functions by clicking the lb, sb, and gf toolbar icons until all of the weighting functions say "unused".

		)									لد - ۱۱	
Exp:1	Seq: PROT	ON	Index: 1	1								
Ér: Ác									Å	J	1.1.	<u>813</u>
						<u>.</u>		**************************************				
TY		**					ž					2
Mouse but	tons: Left - weight	ing Center	- vf/vs_Bigk	t - spectrum	on/off							
vf 4142.591	vs 3057.979	2:lb unused	3:sb unused	4:sbs unused	5:gf unused	6:gfs unused	7:awc 0.034	>			-	₽×

Set the sinebell function to be about double the width of the data in the FID to be kept.



Next click the sinebell shifted icon on the toolbar and shift the sinebell over so that it just covers the FID.



After a good adjustment is found, type **wft** in the command line. This performs a weighted Fourier transform and applies the weighting function. The spectrum may need to re-phased slightly after this.

## Zero-Filling

Zero-filling is handled through the **fn** parameter. The Fourier number, **fn**, is normally set to 'n' or not used. In this case, the number of points (**np**) are the actual data points. To use zero-filling, setting the **fn** to a number larger than the **np** will zero-fill. This number is typically a power of 2. For example type **fn=16k** or **fn=32k** or **fn=np*2** on the command line and the computer will automatically adjust to the correct number. Re-transform (**wft**) the data after changing this parameter to see the result. Using a **fn** smaller than **np** will use fewer than the actual number acquired.

## Measuring Signal-to-noise

Measure the intensity of the largest peak in the spectrum. (Use the cursor to move onto the line and type **nl**.) Then move the cursor over and enclose an area of noise with two cursors. Type **dsn** on the command **line** for a display of signal-to-noise. This command can be found in the Command and Parameter Guide.

For this to be an accurate measure between spectra, use the same peak and the same window of noise each time.

### **Finding Digital Resolution**

Place the cursor near the maximum of the peak to measure. Type **nl dres** on the command line for a display of the digital resolution (width at approx. half height).

### Integration

Display the spectrum of interest on the screen.



Clear any existing integrals by clicking Clear Integrals on the Cursors/Integration Page of the Process tab or by typing cz in the command window.



Click the integral icon on the tool bar to access the integration mode. This will cause a green integration line to appear on the screen.



The different integration icons are shown below.



Click the define integrals icon on the toolbar, then click the left mouse button to cut the green line on either side of the peaks of interest. Clicking the right button will cause the cut to be erased.



To correct leveling and tilt errors in the integral click the level/tilt integration icon



and adjust the integral similar to phasing the spectrum. To set the integral values, put the cursor on one of the peaks and type a numerical value into the Integral Area Box.



Click Set Integral Value and Show integral value to see the integral values on the screen.



### **Baseline Correction**

After integral regions have been specified on a spectrum, the baseline can be corrected to reflect these regions as peaks. To baseline correct, type bc on the command line or click the BC Correct button found on the Display page of the Process Tab. This baseline correction assumes everything in a defined integral region is a peak and flattens the rest. To work properly, everything that is or might be a peak must be in an integral region. To massage this further, see the Varian parameter guide on **bc**. This command can have many modifiers to produce a better correction.



The spectrum below illustrates the use of the baseline correction in a ¹³C 1D spectrum.



# **Peak Picking**

Use the Show/Hide threshold icon to display the threshold level.



Click the Find Peaks button on the Display page of the Process Tab.

Weighting	Display Mode	Axis	Amplitude Scaling	Reference	Baseline Correct	
Display Linear Prediction	<ul> <li>pnased</li> <li>absval</li> <li>power</li> </ul>	<ul> <li>○ Hertz</li> <li>● PPM</li> <li>○ kHz</li> </ul>	<ul> <li>Normalized</li> <li>Absolute</li> </ul>	By Solvent By TMS Cancel	Autofind Integrals BC Correct	
Line Lists Plot Text Output	Screen Posi Full C Left	<b>tion</b> Center Right	Scale Adjust Autoscale + -	By Cursor Reference cursor to 0.00	Find Peaks Peak Threshhold Find Peaks	
	Display Arr Horizontal Vertical C	r <b>ays</b> Label Justom	Display offsetshorizontal-50.3vertical0.0	Nearest Line	Mark at Cursor Clear Marks	

This will display the peak frequencies on all of the peaks which cross the threshold.



Readjust the threshold by clicking the Show/Hide Threshold icon and moving the yellow line with the left mouse button.



Click the Find Peaks button to display the peak frequencies.

### **Basic 1D Plotting**

The simpliest way to print is with a list of commands. These commands are listed below. The following commands should be typed in the vnmr command window to print the following options. They should be typed in a string separated by spaces ending with the command page. All commands are optional except for page. The spectrum will be plotted at a width of 250 mm, unless you specify different. You can change the width of the plot by typing wc=200. This will change it to 200 mm.

ist)
ly

You may also use the buttons on the Process>Plot page. Simply click the options that you would like to print and then click Plot Page.

### Making a PDF Files

To make a pdf file from your spectra, choose the features that you would like to print, then simply click Plot Preview instead of Plot Page.

P					
Basic	Automatic Plot Page	Plot Setup	Manual Plot	FIDs	
Default		As Displayed 🔹	Plot Spectrum	Plot FID	
vveighting	Auto Plot Preview		Dist Constant America	Plot FID Array	
Display Mara 1D		Parameters	Plot Spectrum Array	Plot FID Scale	
Integration	Send this plot to:	Basic (left) 💌	Plot Spectrum Scale		
Cursors/Lino Lists	NMR_room_plot [b+w]	Plot Parameters		Plot Pulse Sequence	
Plot	Hz to mm: 8.4	Integrals	Micc		
Text Output			MISC	Clear Plot	
	Screen Position	Show Plot	Plot Text		
	Full Center	51104	Plot Molecules	Plot Page	
	Left Right	Peak Frequencies		Thot Tage	
	Autoscale	•	Plot Logo		
	WYSIWYG	Show + - Plot	(	Plot Preview	)

Adobe Acrobat will open and show you your pdf file. Click File>Save As> to save the pdf file.

olications Places S	ystern 🝺 Ә 🍭 🍣	ê 💼		
		Vnmrj		
jile <u>E</u> dit ⊻iew E <u>x</u>	periments Acquisition	Process Iools Help		
🎦 🗁 🗔 🚱 🖸	0			
		tmpplot.pdf - Adobe Reader	-0×	
Ele Edt 1	Jew Document Jools	<u>Window</u> Help	×	
	- 4 1 /1	8 8 775N ·   🖶 🚼   Evil - ·		
	5td Proton parameters			
2	Pile: home/vomel/vomesys/de	ta/test.did		
	Pulse Seguence: sipul			
	Solvent: aretone Ambient temperature			
	Operator: vest			
	Mercury-100 *Tarias-ING*			
	Balas. delay 1.000 and			
	Arg. time 0.000 sec			
	Midia 4000.0 Ma Single scan			
	DATA PROCEDUIND			
	Total time 0 min, 14 sec	🐔 Save a Copy 🗙		
		Nama-		1
		Banke. Susan bui	000	S.
		Save in folder: 🙀 vnmr1 🗘		
			8.1	
		P Browse for other folders		
		X Cancel Save	Eu	_
				-
			trum	
			trum A	
		1	trum Se	cale
		AN YMAN		
	فيتعاجب المتحاط والمتحدة	المتحفيدة المجلي المناج المتحاج المتحاج المحاج ا	Positio	on
	de l'alles le stratigies et	feedwikests, tas dealers a feater.	Cer	ner
	8.6 8.5	8.4 8.3 8.2 8.1 8.0 7.9 7.8 7.7 ppm	oscale	
			STWYG	_
0				
G			112	/mm
Temp Spi	n Lock Pro	De la file anno in file deserver interneter et		
<u>111</u> C 20	Hz 37.2 Susan	Prot the saved to the /vmm//mp/tmpptot.ps.		
Vomrt	The transist of	f - Adobe Reader		
	- unppior pr			

### Converting 1D Data to ascii Format

Spectra can be plotted by excel or other programs if it is first converted to ascii format. To convert the FID, first process, phase and reference the data. Take note of the frequency at which it was collected. (Type sfrq? on the command line and write down the number). Type **writexy** on the command line. There will be no indication or message that anything was written, however, the file should have been created inside of the experiment directory. For example, if the vnmr1 group was working inside of exp1, then the output file named xytrace.1 would be found inside of /home/vnmr1/vnmrsys/exp1/ directory. The file will have two columns. The first column is frequency and the second column is intensity. After opening the file in excel, divide the first column by the sfrq value noted earlier to convert it into ppm. Plot data as desired.

### XIII. Finding a 90 Degree Pulse Width

#### Introduction

What is the 90 degree pulse width? The radio frequency pulse is described by its power and duration (time). Before the pulse the proton "spin" is oriented in the z direction and it cannot be observed in the spectrum. The rf pulse "flips" the spin into the xy plane so that it can be observed. The maximum signal is seen when the spin is completely in the xy plane without any z component. The length of time that this takes is called the 90 degree pulse. In the spectra below the length of the pulse is varied from 1 to 40 microseconds. The most accurate way to determine the 90 degree pulse width is to find the 360 and divide it by 4. This value is directly affected by the pulse power. As you increase the power, then the 90 degree pulse will decrease. Normally a power (tpwr) between 54-60 is used.



### Instructions

Collect a ¹D 1H spectrum. Expand the region around one peak as shown in the spectrum below.



Take note of the previously calibrated 90° pw and power level and record them below. These values can be found on the acquisition page of the Acquire tab.

	Start Acquire	Process Show Tin	ne Acquire S	top MoveSW Est. p	w90 Sequence	
	Default H1 Acquisition Pulse Sequence <del>Channels</del> Flags Future Actions Overview	PROTON Dat Spectral width Acquisition time Complex points	6410.3 Hz ▼ 2.556 sec ▼ 16384	Display Sequence Excitation Relaxation delay First pulse Inter-pulse delay Observe Pulse	Arrays 1.000 sec	or 0 degrees
3. .t		Scans Requested Completed Steady-State	8 0 0 v off	Receiver Gain: Calibration: pw90	30 Auto 6.60 at Power	57

Create an array of values for the pulse width (pw). Open the array pop-up window by clicking the arrays button found on the Default 1H tab of the Acquire panel. (Alternatively type array on the command line).

<	Steet Acquire Default H1 Acquisition Pulse Sequence Channels Flags Future Actions	Acquire         Acquire           Process         Acquire           Spectral Width (ppm)         14           Downfield [14.0         Upfie           Pulse Width (degrees)         45           Enter pulse angle         Relaxation Delay (sec)         1           Number of Scans         1           Spin         26         Tune method	ire ê Transform Stop St → -2 ▼ Transfor Id -2.0 ↓ Line Broz 45 ♀ Plotting t ↓ Param ♥ Plotting t Plotting t	ow Time Sequence	Arrays	
	Temp 20.0 C	Spin Lock UHz 80.5 vuto	Probe	ldle	exp2: Acquisition complete	

The array pop-up window is shown below:

_	Application	s Actions	2 6 6 6 6				
$\mathbf{v}$			Array	Paramete	er		
	Param Nan	ne	Description		Size	Order	On/Off
		Array Size:		Total 7	"ime:0:0:1		
			UnArray		New Ar	ray	
	Active Par	am:	None	Curren	it Value:		
					Position		Value
				1	FOSICION	NA	value
	Array Size						
	First Value	2:					
	Increment						
	Last Value						
	Inc. S	ityle	None				
	Rand	omize					
		Edit	Undo Clos	e /	Abandon		

Create the array by doing the following:

Type pw in the Param Name Box Click New Array Fill in the Array Size, First Value, and Increment Boxes Click enter on the keyboard.

The array needs to go past the  $360^{\circ}$  pw. This number can be estimated by taking the previously calibrated  $90^{\circ}$  pw value and multiplying it by 4. You will want your array to go past this point. In this example the  $90^{\circ}$  pw is 11.25 microseconds. The array must go past 11.25*4 (45 µseconds), therefore pw is arrayed from 1 to 55 µsecond in steps of 1 µsecond. After the array is created click the Close button at the bottom of the box.

		Arra	y Paramet	er		
Param Name		Description		Size	Order	On/Off
pw		Pulse width		55	1	On
Ar	ray Size:5	5	Total	Time:0:2:54		
		UnArray		New A	rray	
Active Param	ĸ	pw	Currei	nt Value:	5.95	
				Position	Va	alue
Array Size:		55	1		1 2	
First Value:		1	3		3	-
Increment:		1	5		5	
Last Value:		55	6 7		6 7	
Inc. Styl	e	Linear	8		8	
Random	nize		10		10	
			11		11	
			13		13	

An appropriate delay, number of scans, absolute intensity must be now set and the experiment started. A delay of 2 seconds with 1 scan is chosen in the example. These are set by typing:

	d1=2 nt=1 ai ga on the comman	d line.
	📥 Applications Actions 🗾 🤗 🥸 🖏 🐻	Tue Nov 20, 09
	Vnm Vnm	J = = ×
	<u>File Edit View Experiments Acquisition Process</u> Tools Help	
	: 📩 📂 🗔 🚱 🔘 🔘	- 🖸 <b>- IL</b> - M
		Q. J. <u>M. 11</u> <del>M.</del> <b>Z</b> O ×
	exp2: Acquisition complete	
/		<b>_</b>
(	✓ d1=2 nt=1 ai ga	
	Exp:2 Seq: Proton Index: 1	

When the experiment is completed display one spectrum from the array and phase it. In this example the third spectrum is displayed and phased. This is done by typing **ds(3) full aph** on the command line.



The 90° pw is calculated by taking the 360° pw and dividing it by 4. The 360° pw is identified as the second null point in the curve. The numbers can be displayed under each spectrum by typing dssl (not shown). In this example, the 45th spectra is the 360° pw. This corresponds to 45  $\mu$ seconds.



To find a more accurate  $90^{\circ}$  pw create a second array that spans 4  $\mu$ seconds around the  $360^{\circ}$  pw with 0.25  $\mu$ seconds increments. In this example, the second array would be from 43 - 47  $\mu$ seconds with an increment of .25  $\mu$ seconds.

## XIV. Basic ¹³C 1-D NMR

# Proton Decoupled ¹³C Experiment: ¹³C{¹H}

Experimental Set-up:

VNMRS-400 and VNMRS-500: Be sure that the probe is tuned to 1H on channel 1 and 13C on channel 2. Mercury VX-300: the probe is already tuned.

If the proton signals of your compound are either all aromatic or all aliphatic, the default parameters for 13C 1-D NMR may not give good enough signal to noise. In this case, the instructions below can be followed to improve your data.

First acquire a proton spectrum of the sample in exp1. Move the transmitter offset to the center of the proton signals by placing the cursor in the middle of the peaks and typing **movetof** in the command window. This value may then be found by typing tof? in the command window.





Note the **tof** (transmitter offset). This will become the set value of dof in the carbon experiment.

Join experiment 2 by typing jexp2 in the command window. If experiment 2 doesn't exist, create it by typing cexp(2). If experiment 2 is locked, unlock it by typing unlock(2). Set up the carbon experiment by choosing the carbon experiment under the Experiments drop down menu.



Set dof to the center of the proton peaks (tof) that was determined above. Type dg and the parameters will be shown in the text output (or overview) window of the process tab.

dof=996.6 dg									
<b>_</b>									
Exp:2 Seq: C	ARBON	Index:	1						
Start Acquir Pr	ocess	Transform	Autoprocess	Display	/ Spectrum	Clear Scre	en <b>e</b>	Cancel	ф.;
Start Acquir Pr Basic	ocess	Transform .	Autoprocess	Display	/Spectrum SPE	Clear Scre	een PF	Cancel	<b>.</b>
Start Acquir Pr Basic Default	ocess seqfi	Transform QUISITION 1 s2pu1	Autoprocess TRANSM: tn	Display LTTER C13	/Spectrum SPE temp	Clear Scre CIAL 25.0	en Pf 1b	Cancel ROCESSING 0.50	. <del>П</del>
Start Acquir Pr Basic Default Weighting	ocess seqfi sw	Transform QUISITION 1 S2pu1 25510.2	Autoprocess TRANSM: tn sfrq	Display CITTER 100.582	/ Spectrum SPE temp spin	Clear Scre CIAL 25.0 not used	een Pr 1b sb	Cancel ROCESSING 0.50 not used	<del>р.</del>
Start Acquir Pr Basic Default Weighting Display	ocess seqfi sw at	Transform QUISITION 1 S2pul 25510.2 1.285 55536	Autoprocess TRANSM: tn sfrq tof town	Display C13 100.582 1531.4	/Spectrum SPE temp spin gain bet	Clear Scre CIAL 25.0 not used not used	een Pr 1b sb gf 3wC	Cancel ROCESSING 0.50 not user not user	
Start Acquir Pr Basic Default Weighting Display More 1D	ocess seqfi sw at np fb	Transform QUISITION 1 \$2pu1 25510.2 1.285 65536 17000	Autoprocess TRANSM: tn sfrq tof tof tpwr pw	Display ITTER 100.582 1531.4 53 3.650	/Spectrum SPE temp spin gain hst pw90	Clear Scre ZIAL 25.0 not used not used 0.008 7.300	een PF 1b sb gf awc 1sfid	Cancel ROCESSING 0.50 not user not user not user	
Start Acquir Pr Basic Default Weighting Display More 1D Integration	ocess seqfi sw at np fb bs	Transform           QUISITION           1         \$2pul           25510.2         1.285           65536         17000           8         8	Autoprocess TRANSM: tn sfrq tof tpwr pw DECOUR	Display C13 100.582 1531.4 53 3.650 21 FR	/Spectrum SPE temp spin gain hst pw90 alfa	Clear Stree CIAL 25.0 not used 0.008 7.300 10.000	een Pr 1b sb gf awc 1sfid fn	Cancel 0.50 not user not user not user not user	
Start Acquire Pr Basic Default Weighting Display More 1D Integration Cursors/Line Liste	ocess seqfi sw at np fb bs d1	Transform QUISITION 1 \$2pu1 25510.2 1.285 65536 17000 8 1.00	Autoprocess TRANSM: sfrq tof tpwr pw become dn	Display C13 100.582 1531.4 53 3.650 PLER HI	/Spectrum SPE temp spin gain hst pw90 alfa PRESAT	Clear Stree CIAL 25.0 not used 0.008 7.300 10.000 JRATION	PF 1b sb gf awc 1sfid fn	Cancel ROCESSING 0.50 not user not user not user not user FLAGS	
Start Acquir Pr Basic Default Weighting Display More 1D Integration Cursors/Line Lists	ocess seqfi sw at np fb bs s d1 nt	Transform QUISITION 1 \$2510.2 1.225 65536 17000 1.000 1.000	Autoprocess TRANSM: tn sfrq tof tpwr pw DEcound dn dof	Display CITTER (13 100.582 1531.4 53 3.650 01EP H1 996.6	/Spectrum SPE temp spin gain hst pw90 alfa PRESAT atmode	Clear Scre CIAL 25.0 not used not used 0.008 7.300 10.000 JRATION n	Peen PF 1b sb gf awc 1sfid fn il	Cancel ROCESSING 0.50 not user not user not user not user FLAGS	
Start Acquir Pr Basic Default Weighting Display More 1D Integration Cursors/Line Lists Prot	occess seqfi sw at np bs d1 nt ct	Transform           QUISITION           1         \$2pul           25510.2         1.285           65536         17000           8         1.000           000000000000000000000000000000000000	Autoprocess TRANSM: tn sfrq tof tpwr pw pecould dof da	Display CTTER (13 100.582 1531.4 53 3.650 PLER H1 996.6 996.6	/Spectrum SPE temp gain hst pw80 alfa alfa PRESAT wet	Clear Scre 25.0 not used 0.008 7.300 10.000 JRATION n	PP 1b sb gf awc 1sfid fn il in	Cancel O.50 not used not used not used not used FLAGS	,

To determine the number of scans for your available time, type time(hours,minutes). For example, if I have 1 hour and 15 minutes, type time(1,15) dg. This will set the number of scans (nt) to the appropriate number. Type bs=16 ga wbs('wft') in the command window. This will start the experiment and show you the carbon experiment at the completion of every 16 scans. If you would like to stop the experiment before it is complete, click the Stop button or type aa. **Save the data**, send it to the data station and print it.

### Important Parameters:

The parameters dpwr, and dmf are essential for good results. These parameters are calibrated regularly and correct values should load into the parameter set when the experiment is set up. These parameter values along with their calibration dates can be found in the red calibration notebook located next to each spectrometer. The decoupler power (**dpwr**) should stay at or below 40db for safety. The decoupler modulation mode (**dmm**) is the mode the decoupler is using. To decouple protons, use waltz decoupling (**dmm='w'**). When using
waltz or garp decoupling, the decoupler modulation frequency (**dmf**) is set to 4X  $\gamma$ H2. The  $\gamma$ H2 value is regularly calibrated and updated in the probe file.

The spectrum below was collected with 32 scans. The decoupler was on (dm='yyy'). The spectrum shows three singlets at 124.0, 126.7 and 127.5 ppm.



# Integrating the Carbon Spectra

In a proton decoupled ¹³C spectrum all the coupling information between the 1H and 13C has been removed and all the carbons peaks should be singlets. Integration of this spectrum may not be accurate due to the carbon-proton NOE enhancement. The default setting in the carbon pulse sequence accessible by the drop-down menu for the decoupler is "on" for the entire length of the experiment. This is designated by dm='yyy'. Use **dps** (display pulse sequence) to see when the decoupler is off and on. Decoupling the protons in the carbon experiment significantly increases signal intensity, but removes coupling information and the possibility of integrating the spectrum.

# Proton Coupled ¹³C Experiment

Carbon-13 1D NMR can also collected with the C-H coupling observed. This is referred to as a proton-coupled carbon-13 experiment: Different decoupling settings will allow for better or worse signal-to-noise and determines whether the spectrum can be accurately integrated and/or whether the coupling information (splitting) is retained.

<b>Decoupling Mode</b>	S/N	Integration	Splitting
'nnn'	poor	yes	yes
'nny'	fair	yes	no
'yyn'	fair	no	yes

|--|

Set the decoupler to the desired value. The spectrum below was collected with the decoupler off (**dm='nnn'**). This change in decoupler settings resulted in a large decrease in signal-to-noise but contains the coupling information. Many transients were needed to see carbon peaks (1798 scans). Carbon peaks from carbons attached to protons are split. This spectrum can be integrated. The spectrum shows three doublet of doubles at 124.0, 126.7 and 127.5 ppm (J₁=163 Hz, J₂=9 Hz).



# XV. 2D- Experiments

# 2-D Parameter Sets

Typical operation usually collects the normal 1-D in exp1 and 2-D experiments in exp2 or higher. Virtually all 2-D experiments will require a 90 deg. pulse width calibration for setting up the experiment. Also move the spectral window as necessary to observe only the area of interest and set up any water suppression needed.

For homonuclear 2-D experiments like COSY and NOESY, the normal 1-D parameters control the f2 dimension (direct dimension) while the f1 axis (indirect dimension) is controlled by parameters identified with a 1 (i.e. **sw1, lb1, fn1**). In homonuclear experiments, be sure **sw=sw1**. One other very important parameter in 2-D experiments is the number of increments (**ni**). The **ni** specifies the number of "points" in the indirect dimension. Typical values for ni are 128, 256, or 512. Steady state scans are also important in long experiments so be sure to set **ss**. The suggested value is **ss=32**. A manual page (**man**) for the experiment will also point out parameters that need to be set. To access these type the following on the command line:

man('pulsesequence name') or printon man('pulsesequence name') printoff i.e. man('NOESY') To determine the length of time the experiment will run, type **time** on the command line. Be sure the experiment fits into the time block available. To change the length of time the experiment will run, options include increasing or decreasing **ni**, **nt** or **d1** (or **satdly** in a presat experiment). Check the **at** (acquisition time) to check that a typical number of points is being taken. The typical number of points (**np**) taken is 2048 or 4096. These correspond to **at** times of less than 0.5 seconds. The **at** time in some 2-D experiments *must* be kept short due to high power decoupling during this time period.

### General Considerations for doing 2D Experiments (Varian inc, "Two-Dimensional NMR 5 days Series").

1. All 2D experiments will give better results (less t1 noise) if they are run with the spinner turned off. FOR HMQC, HMBC, INADEQUATE and ALL water suppression of gradient experiments it is mandatory that the spinner be turned off.

2. All 2D experiments will give better results if sample temperature is regulated during data acquisition.

3. All 2D experiments will give better results if lock conditions are optimized. This involves setting Z0 more or less exactly on resonance for the lock solvent, adjusting lock power to just below saturation (varies from solvent to solvent), optimizing lock phase and setting lock gain as low as possible (after shimming).

4. The S/N of individual F2 spectra (before second FT) needs to be only about 5:1 in many cases. In partitioning total experiment time between scans per increment (nt) and number of increments (ni), it is probably better to emphasize S/N of individual increments a little more heavily, bearing in mind the square root dependence of S/N on number of scans. Increasing the number of increments (ni) excessively may make linear prediction less advantageous, reducing overall efficiency.

5. It is always advantageous to reduce the spectra width in both dimensions to the minimum practical value, keeping in mind how peaks may fold in.

6. Zero-filling in F1 is universally done. Linear prediction is also quite common, typically being used to extend F1 interferograms by a factor of 2 in length. (If linear prediction is used, the added points must be included in determining how many zero-filling points to add. The F1 window function must also be adjusted to include the added points).

**COSY**-(Correlation Spectroscopy) spectra are 2D proton homonuclear spectra that correlate one-bond J couplings. A basic COSY sequence is comprised of two 90-degree pulses.

# Versions of COSY

gCOSY – gradient version of the upper case COSY, faster if abundant sample, allows down to nt=1. With nt=1 and ni=128, a gCOSY can be run in as little as three minutes. Gradient COSY setup is detailed below. In gCOSY, phase=1. No presaturation option available.

*cosy* (lower case) – older version of COSY, no gradients, sets up relay cosy (relayh)

Set up by typing **cosy** on the command line and setting **pw**, **sw=sw1**, **nt=4** (min.), **d1**, and check >**time**. **relay=0** for absolute value cosy. (set **relay** and **tau** for relay or long-range COSY)

*cosyps* (lower case) – phase-sensitive COSY, no gradients, similar to **cosy**, but allows solvent presaturation option and data is phase sensitive. Minimum **nt=4**.

#### DQF-COSY

The double-quantum filtered COSY is used to examine coupling constants in detail. In general, this experiment takes much longer than a regular COSY but can provide much more information when measuring coupling constants in multiplets.

# gCOSY—Experimental Set-up

Put your sample into the magnet. Lock, shim, tune and collect a proton experiment in experiment 1. Calibrate a 90 degree pulse width.

Join experiment 2 by typing jexp2 in the vnmrj command window. Collect a proton spectrum and adjust the sweep width, gain, tof, and set the pw90 to the calibrated 90 degree pulse width.

	Start Acquire I	Process Show Tin	ne Acquire	Stop MoveSW	Est. pw90	Sequence	
(	Default H1 Acquisition	PROTON		Display Sequence	Arr	rays	
	Pulse Sequence	Data		Excitation			
	Channels	Spectral width	6410.3 Hz	<ul> <li>Relaxation</li> </ul>	delay 1.000	sec 💌	
	Flags	Acquisition time	2.556 sec	<ul> <li>First pulse</li> </ul>	0.00	us 🔻 or O	degrees
	Future Actions	Complex points	16384	Inter-pulse	delay 0.000	sec 💌	
	Overview		1	Observe Pi	ulse 3.30	us 🔻 or 45	degrees
_		Scans		00000000	inge large		degrees.
۶.		Requested	8	D	20		
_		Completed	0	Seceiver Ga	n. po		
ıt		Steady-State	0 🔽 off	Calibration:	pw90 6.60	at Power 57	/
-							

Change the experiment into a gradient COSY experiment by using the Experiment drop down menu at the top of the page.



Click the gradient COSY button. Notice that the sequence at the top now reads gCOSY.



Look at the parameter list on the Overview page of the Acquire tab. Type dg to refresh this parameter list. Make sure that pw, tpwr, sw and tof are the same as previously noted. If not, change them by typing the abbreviation and value (i.e. sw=4803). The sweepwidth for the second dimension should equal the sweepwidth for the first dimension in this experiment. Therefore sw1 should equal sw.

			А		l				В	
	Start Acquire F	rocess	Show Time	Acqu	uire Stop				Sequence	
	Defaults	ACQ	UISITION	TRA	NSMITTER	GRAD	IENTS	P	ROCESSING	TI
	Acquisition	seqfil	gC0SY	tn	H1	gzlvlE	5102	sb	-0.075	lext
	Pulse Sequence	SW	6410.3	sfrq	399.965	gtE	0.001000	sbs	not used	D I I I I I I I I I I I I I I I I I I I
	PRESAT	at	0.150	tof	399.9	EDratio	1.000	fn	2048	Basic
	WET	np	1924	tpwr	57	gstab	0.000500	ZD	PROCESSING	
	VVL I	SS	32	pw	6.600	hsglvl	6120	sbl	-0.020	Array
	Channels	dl	1.000	DE	COUPLER	hsgt	0.002000	sbsl	not used	
	Flags	dZ	0.020	dn	C13	FL	AGS	proc1	1p	Channels
	Future Actions	nt	1	cim	nnn	hs	nn	fn1	2048	
2	Overview	ct	0	S	PECIAL	sspul	У			Solv.Supp
- I	Overview	ZD A	CQUISITION	temp	not used	SA	MPLE			
		sw1	6410.3	spin	0	date Mar	ch 8, 1993			Shime
		ni	128	gain	30	solvent	dZo			511113
Jt.		PRES	ATURATION			sample				Chan
		satmode	e n							Crear
		wet	n							
_		Detai	ls in dass							

Type nt=1 ss=4 ni=1 ga and watch the remote status box for a ADC overflow error. Raise the gain value until the signal overflows, then set the gain to 5 less than the overflow value. (The maximum value is 60).

Determine the necessary number of transients (nt) by observing how many transients are required to see the proton signal.

Type nt = required number as determined above, ss=32 ni=128 Type time to see how long the experiment will be. Click the green Acquire button.

Save the experiment when it is completed.

**NOESY** (Nuclear Overhauser Effect Spectroscopy) spectra provide information about protons that are 5 Angstroms or less apart in space. The information is through space and not through bond, like a COSY. The presence of a NOE peak is direct evidence that two protons are within 5 Angstroms through space. The absence of a NOE peak between protons does not necessarily mean that they are not within 5 Angstroms since other factors can reduce a NOE peak even if the protons are close in space. A mid-size molecule (~1000-1500MW range) may have NOEs that are close to zero and a ROESY may be required to see them. Large molecules generally give better NOEs at higher field, but small molecules may actually give better NOEs at lower field. A 2-D NOESY of a small molecule will have cross peaks of opposite phase to the diagonal. A 2-D NOESY of a large molecule will have cross peaks of the same phase as the diagonal. Theoretically, these experiments should be symmetrical, but it is typical to see more intense peaks on one side of the diagonal than the other.

#### Experimental Set-up

Put your sample into the magnet. Lock, shim, tune and collect a proton experiment. Calibrate a 90 degree pulse width; adjust the sweep width and tof. Convert the proton data experiment into a NOESY experiment using the drop down menu.



Check that **pw** is set to the calibrated 90-degree pulse width, set **d1** and **np**. Set the **sw1=sw**. Set **phase=1,2**. Type nt=1 ss=4 ni=1 ga and watch the remote status box for a ADC overflow error. Raise the gain value until the signal overflows, then set the gain to 5 less than the overflow value. (The maximum value is 60).

Determine the necessary number of transients (nt) by observing how many transients are required to see the proton signal. The minimum nt in this experiment is 2.

One parameter requiring forethought is the **mix** time. The **mix** is usually determined by the size of the molecule under study. Small molecules require longer mix times, 0.5 to 0.8 seconds (suggested **mix=0.5**). Large molecules generally range from 0.05-0.3 seconds (suggested **mix=0.15**).

Set the **ni and ss=32** then type **time** to determine experimental time. Adjust the delays, transients and/or increments to fit the experiment into the time available. Click the acquire button to start the experiment.



To do simple processing of these spectra type **setLP1 gaussian wft2da** on the command line. Further phase correction and manually processing may be required to improve the spectral quality. These corrections are described in the 2D processing and printing section of this manual.

#### Indirect Detection Experiments

Indirect detection is the detection of a heteronucleus through direct detection on proton signals. Detection on proton allows for a much higher sensitivity. Keep in mind that the amount of sample required for a heteronuclear experiment is still much larger than a purely proton experiment.

**HSQC** (Heteronuclear Single Quantum Coherence) A proton-carbon HSQC will detect all carbons with a proton attached. Carbons with no proton attached will not appear in the 2-D spectrum.

#### Experimental Set-Up

Put your sample into the magnet. Lock, shim, tune (to both proton and carbon) and collect a proton experiment. Calibrate a 90 degree pulse width; adjust the sweep width (sw) and tof. Take note of the parameters pw90, tpwr, sw and tof and check that they are carried over into the HSQC experiment. Convert the proton data experiment into a HSQC experiment using the drop down menu.



#### Parameters:

Set **sw1**=carbon sweep width. Use a carbon sweep width large enough to encompass all the carbons attached to protons (don't worry about quaternary shifts). Set the **dof**=decoupler offset for carbon. If certain what **dof** to use, a chart of transmitter offsets for carbon is found in the red calibration notebooks located next to the spectrometer's computer. **pwx** is the 90 degree pulse width of carbon at the power level **pwxlvl**. This is normally calibrated with a standard sample. The value is found in the probe file. The decoupling modulation frequency, **dmf**, should be equal to the necessary window for decoupling. **dmf** is calibrated with a standard sample and is calibrated at power level **dpwr**. **dmf** is equal to 1/(90 degree pulse width). These values are found in the probe file. Set the **at**= acquisition time. Make this time short (<0.2 seconds) due to high decoupling power on carbon during this time. Set **phase=1,2** for HSQC. The one-bond coupling is selected by **j1xh**. Type nt=1 ss=4 ni=1 ga and watch the remote status box for a ADC overflow error. Raise the gain value until the signal overflows, then set the gain to 5 less than the overflow value. (The maximum value is 60).Set the **ni and ss=32** then type **time** to determine experimental time. Adjust the delays, transients and/or increments to fit the experiment into the time available. Click the acquire button to start the experiment.



To do simple processing of these spectra type **setLP1 gaussian wft2da** on the command line. Further phase correction and manually processing may be required to improve the spectral quality. These corrections are described in the 2D processing and printing section of this manual.

# HMBC (Heteronuclear Multiple Bond Coherence)

HMBC is a heteronuclear 2-D experiment that will pick up carbons (non-isolated) without a proton attached. This experiment can be very useful to see quaternary peaks when there isn't enough sample for a 1-D carbon. The HMBC set up is very similar to HSQC except there is no decoupling on carbon.

# Experimental Set-Up

Put your sample into the magnet. Lock, shim, tune (to both proton and carbon) and collect a proton experiment. Calibrate a 90 degree pulse width; adjust the sweep width (sw) and tof. Take note of the parameters pw90, tpwr, sw and tof and check that they are carried over into the gHMBC experiment. Convert the proton data experiment into a gHMBC experiment using the drop down menu.



#### Parameters:

Set **sw1**=carbon sweep width. Use a carbon sweep width large enough to encompass *all* possible carbons in the molecule. Set the **dof**=decoupler offset for carbon. If certain what **dof** to use, a chart of transmitter offsets for carbon is found in the red calibration notebooks located next to the spectrometer's computer. **pwx** is the 90 degree pulse width of carbon at the power level **pwxlvl**. This is normally calibrated with a standard sample. The value is found in the probe file. Set the **at**= acquisition time. Acquisition time can be longer here since the decoupler is off during this experiment. Set **phase=1,2** for gHMBC. The one-bond coupling is selected by **j1min and j1max** and the multiple bond by **jnxh**. Type nt=1 ss=4 ni=1 ga and watch the remote status box for a ADC overflow error. Raise the gain value until the signal overflows, then set the gain to 5 less than the overflow value. (The maximum value is 60).Set the **ni and ss=32** then type **time** to determine experimental time. Adjust the delays, transients and/or increments to fit the experiment into the time available. Click the acquire button to start the experiment.

Start Acquire F	Process	Show Time	Acquire	e Stop				Sequence
Defaults	ACOUT	STTTON	TRANS	MITTER	. н	MBC	PRO	DCESSING
Acquisition	SW	6410.3	tn	H1	j1min	130.0	sb	-0.075
Pulse Seguence	at	0.150	sfrq	399.965	j1max	165.0	sbs	not used
DDECAT	np	1924	tof	399.9	jnxh	8.0	fn	2048
PRESAT	bs	32	tpwr	61	GRAD	IENTS	2D PF	ROCESSING
WEI	SS	32	pw	8.800	gzlv1E	4524	gf1	0.015
Parameters	d1	1.000	DECO	UPLER	gtE	0.002000	gfs1	not used
Channels	nt	4	dn	C13	EDratio	3.977	proc1	1p
Flags	ct	0	dof	1028.5	gstab	0.000500	fn1	4096
Futuro Actione	ZD ACQ	UISITION	pwx1v1	53	hsglvl	5424	2	SAMPLE
Future Actions	sw1	24140.0	pwx	7.600	hsgt	0.002000	date Ma	arch 8, 1993
Overview	ni	200	dm	nnn	SPE	CIAL	solvent	t D20
	phase	arrayed	decwave	9	temp	not used	sample	
	PRESAT	URATION	dmf	29412	spin	0		
	satmode	n	dpwr	31	gain	30		
	wet	n			pw90	8.800		
	Details	in dgss			sspul	У		

To do simple processing of these spectra type **setLP1 gaussian wft2da** on the command line. Further phase correction and manually processing may be

required to improve the spectral quality. These corrections are described in the 2D processing and printing section of this manual.

XVI. 2D Processing and Printing

# 2D Tool Bar

Icon	Description
	Box/Cursor
6	Show Full spectra
	Zoom In Zoom Out Zoom Mode
~~ <u>~</u>	Pan/Stretch Mode
~	Trace
L.L.	Scale
é	Projections
3	Redraw
2	Rotate
	Raise vertical scale by 20% Lower vertical scale by 20%
14	Phase Mode
	Peak Picking
3	Return

# Manipulating the 2D data

Click on the appropriate icon on the tool bar to expand, zoom, and change the vertical scale of the spectra. Type dconi on the command line to display the peaks using the color map in the interactive mode. Type dpcon to display the peaks in contours. To show more contours than is automatically displayed, use a command that specifies the number of contours and spacing. To display this, **dpcon(20,1.3)** will display 20 contours at a spacing of 1.3. Other variations of this include **dpcon('pos',20,1.2)**. The modifiers 'pos' and 'neg' will display only the positive or only the negative peaks in a 2-D spectrum. To make the contour display interactive, use a variation of **dconi** such as **dconi('dpcon','pos',20,1.3)**. To plot more contours, **pcon(20,1.3) pltext page**. Adjust the number of contours and spacing for the data. The 'pos' and 'neg' modifiers can also be used in the plotting.

#### Processing 2D Data Sets Manually

This routine applies to most phase-sensitive 2-D data sets. COSY requires a slightly different approach and often uses a pure sinebell function.

Start by processing the first increment with by typing **wft(1)** Phase this spectrum as any 1-D.Then add a weighting function to this data by typing **wti** on the command line enabling the use of the interactive weighting screen to add a weighting function. The most commonly used function at this step is a gaussian function **[gf]**.



Then process the data with by typing **wft1da**. After the data is processed, the f1 traces will be on the screen. Select a trace with the cursor. Try to use something other than solvent. Then weight this trace with **wti**.



This weighting function can also be a gaussian function **[gf]**, but often a shifted sinebell is used. To put in a shifted sinebell, start with **[sb]** and move the cursor to produce a sinebell curve that is about twice the width of the interferogram.

Then select **[sbs]**. Using the cursor, shift the sinebell back so the maximum starts at the left side of the interferogram. Process the 2-D data with >wft2da



#### Manual Phase Correction

Phase errors in phase-sensitive 2-D data sets can often be seen near the diagonal where the peaks may be streaked positive and negative. Display the entire 2-D spectrum. Choose 3 traces containing cross peaks near the top, middle and bottom of the spectrum and note the index #'s of each trace. The index # of the trace can be seen in the top window of VNMR next to the seq. and exp #.





Type **r1=***index***# r2=***index***# r3=***index***#** on the command line. The values r1-r3 are place holders in VNMR. Then display the first trace: Type **ds(r1) and** Phase

the 1-D trace. Display the  $3^{rd}$  trace. **ds(r3)** Click  $\checkmark$ . Click the left mouse button on both sides of the spectrum to accept the previous phase changes. Then phase this trace. Go back to the first trace by typing **ds(r1)** Continue phasing r1 and r3, clicking to accept phase changes in-between until they are both phased. Look at the middle trace to check by typing **ds(r2)** Then go back to the 2-D

spectrum by typing **dconi.** If necessary, rotate the axis by clicking the *icon* and repeat the phase correction procedure.

#### **Referencing the 2D Data**

Place the cursor on the contour which you wish to use as your reference. If the spectra is homonuclear, then type rl(#p) rl1(#p) dconi. If the spectra is heteronuclear then type rl(#p) rl1(#d) dconi. Alternatively, the spectra can be referenced by filling in the reference boxes on the Default page of the Process Tab.

Start Acquire Pro	ocess Transform Autoprocess	Display Spectrum Clear Screen	Cancel
Basic Default Weighting Display More 2D Integration Cursors/Line Lists Plot Text Output	Transform         FT 1D - 1st Increment         Transform F2         Full 2D Transform         FT Data Size         Acq Pts         ✓ F1         Ik         ✓ F2         512         300         Transform Coefficients         100000-10         Weighting:         F1       gaussian         F2       gaussian         Save Current Process/Dis	Display         Display ZD         Display Trace         Projections       Full Screen         AutoScale 2D         Trace       ● F1 ● F2         Axis       Display Mode         F1       Phased         F2       PPM         Unear Prediction       ✓         ✓       F1         Auto LP       F1         ↓ F2       Auto LP         splay parameters       fin FID directory/	Display 1D # 424 Display Text BC Correct (F1, F2) DC Correct (F1, F2) Reference F1 by Solvent Paterence F1 by Solvent Set F1 cursor to: 0.00 Hz Set F2 cursor to: 0.00 Hz Set F2 cursor to: 0.00 Plat

Place the cursor of the desired peak, set the units to ppm, enter the correct value into the boxes, enter on the keyboard, then redraw the spectra to observe the referenced spectra.

Basic Default Weighting Display More 2D Integration Cursors/Line Lists Plot Text Output	Transform         FT 1D - 1st Increment         Transform F2         Full 2D Transform         FT Data Size       Acq Pts         V F1       128         V F2       512         Transform Coefficients       100000-10         Weighting:       F1         F2       ✓         Save Current Process/Different	Display         Display 2D         Display Trace         Projections         Full Screen         AutoScale 2D         Trace       ● F1 ● F2         Axis       Display Mode         F1       PF1 ● F1 ● F2         Axis       Display Mode         F1       PPM ▼         Phased ▼          Linear Prediction       ▼         ♥ F1       Auto LP F1         ■ F2       Auto LP F2         isplay parameters       In FID director	Display 1D # 147 Display Text BC Correct (F1, F2) DC Correct (F1, F2) DC Correct (F1, F2) Reference F1 by Solvent Reference F1 by Solvent Reference F1 by Solvent Set F1 cursor to: 39.51 PPM Set F2 cursor to: 2.50 PPM Flot
-----------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# Plotting the 2D Data

The data can be easily plotted using the buttons on the Process >Basic page or the Process>Plot page. Choosing Plot Preview will allow the option of saving the file as a pdf, while choosing Plot Page will send the printout to the printer.

Start Acquire Pro	cess Transform Autoprocess	Display Spectrum Clear Screen	Cancel
Basic Default Weighting Display More 2D Integration Cursors/Line Lists Plot Text Output	Sample: Solvent: dmso sample owner: vnmr1 Comments: STANDARD PROTON PARAMETERS Display Spectrum	Process Options FT data size: ▼ ✓ F1 linear pred. 4*ni ▼ Weighting (F2): gaussian ▼ Weighting (F1): gaussian ▼ Process Save Current Proce	Plot Options arameters: Basic • Plot 1D (F2): spectrum • (F1): spectrum • Plot Plot Plot Preview ess/Display parameters

Start Acquire Pro	cess Transform Au	itoprocess Displa	ay Spectrum Clear	Screen Cancel	
Basic	Automatic Plot	Contour Plot		Plot Pulse Sequence	1D Spectr show plot
Default	Auto Plot Preview	Positive	show plot	Parameters	Plot Spectrum Scale
Display		Negative	show plot	Basic Parameters 🔹 🔻	FID show plot
More 2D	Send this plot to	Both	show plot	Plot Parameters	Plot FID Scale
Integration	Plotter Type: B+W	Stacked Plot	show plot	Misc	The search
Eursors/Line Lists				Plot Text	Clear Plot
Text Output	Screen Position	Projection	Horiz. plot	Top Spectrum	
	Lett +x -x		vert. plot	Side Spectrum	Plot Page
	Right +y -y	Trace	show plot	Plot Molecules	Thorage
	Center "X /X	Trace Axis	● F1 ○ F2	Plot Logo	
	Projections	Hz to mm:	6.0	Plot 2D Frequency List	Plot Preview

#### A.1. 300 NMR Practical Assignment #1 Collecting a 1D proton NMR spectrum

Student Name: ______

1. Collect a ¹D 1H spectrum of 2-Ethyl-1-indanone. Follow the instructions in Section VIII of this manual. Use all default parameters.

2. Reference the chloroform peak to 7.27 ppm. Integrate each peak and set the peak at 7.76 ppm to 1 proton. Print the spectrum from -1 ppm to 10 ppm with the integration values at the bottom of the page. Print a second spectrum from -1 ppm to 10 ppm with the peak frequencies listed on the page.

3. Save this spectrum into your personal directory on the 300 NMR computer. File Name:

4. Transfer this file to Data Station 1. Make sure that it is put inside of your personal directory of Data Station 1. The instructions for transferring data are found within Section V of this manual.

5. Collect a ¹D 1H spectrum of 2-Ethyl-1-indanone. Use a delay (d1) of 2 seconds, an acquisition time (at) of 4 seconds, and a collection window (sweepwidth) of 0 - 9 ppm and 32 scans. The instructions for changing these parameters are found in section XI of this manual.

6. Reference the chloroform peak to 7.27 ppm. Integrate each peak and set the peak at 7.76 ppm to 1 proton. Print the full spectrum with the integration values at the bottom of the page. Print a second spectrum with the peak frequencies listed on the page.

7. Save this spectrum into your personal directory on the 300 NMR computer. File Name: ______.

8. Transfer this file to Data Station 2. Make sure that it is put inside of your personal directory of Data Station 1.

9. Make a copy of this page, attach the four requested printouts to it and hand it into the NMR facility staff.

### A.2. 300 NMR Practical Assignment #2 Collecting a 1D ¹³C{1H} NMR spectrum

Student Name:

1. Collect a proton decoupled ¹³C spectrum of 2-Ethyl-1-indanone following the instructions below.

a. Collect a ¹D 1H spectrum of 2-Ethyl-1-indanone. Follow the instructions in Section VI of this manual. Use all default parameters.

b. Change the 1H experiment into a carbon experiment by clicking the carbon experiment in the Experiments drop down menu.



c. Set a large number of scans by typing nt=100000 bs=16.

						v	nmrJ		
<u>F</u> ile	<u>E</u> dit	⊻iew	E <u>x</u> periments	<u>A</u> cquisition	Automation	<u>P</u> rocess	Tools	<u>H</u> elp	
*	1	- 📀	$\bigcirc$ $\bigcirc$						
:		2							
	nt=10	0000 k	os=16						

Start the experiment by clicking the green acquire button.

ile <u>E</u> dit <u>V</u> iew	E <u>x</u> periments <u>A</u> co	quisition Automation	<u>P</u> rocess <u>T</u> ools <u>H</u> elp			
티 📂 🗖 😨						W - JL -
	br=16					
1111-100000	03=10					
Exp:3 Sec	: CARBON I	ndex: 1				
					2.9	13
	1.0					۸ ۵
	1.0					•\/\∕~
)ec ////////////////////////////////////						
	A			8		c
Start Acquire	Process Show T		MoveSW	Secuen		¢
ocare respanse	The second second			Sequen		
Default C13	CARBON		Display Sequence	Arrays		
Pulse Seguenc	e Data		Excitation			
Channels	Spectral width	25510.2 Hz 🔻	Relaxation delay	1.000 sec 🔻	-1	
Flags	Acquisition time	1.285 sec 💌	First pulse	0.00 us 💌	or 0 degrees	
Future Actions	Complex points	32768	Inter-pulse delay	0.000 sec 💌		
Overview	Scone		Observe Pulse	2.85 us 💌	or 45 degrees	
	Requested	256				
	Completed	0	Receiver Gain:	30 🗹 Auto		
	Steady-State	0 P off	Calibration: pw90	) 5.70 at Power	56	

The data may be processed at the completion of every block of data. This was defined to be 16 scans when bs=16 was typed in the command line above. Type wft on the command line to process the data so that the carbon peaks can be observed. Typing the command wbs('wft') on the command line should cause the spectra to process automatically at the end of each block. Click the red Stop button when the signal to noise is satisfactory.

2. Reference the chloroform peak to 77.23 ppm. (See Section XII of this manual)

3. Save this spectrum into your personal directory on the 300 NMR computer. File Name:

4. Make a pdf file of the spectrum with the peak frequencies printed on the spectrum. Save the pdf file in your personal directory on the 300 NMR computer. (See Section XII of this manual)

PDF File Name: ______.

5. Make a copy of this page and hand it into the NMR facility staff.

# A.3. 400 NMR Practical Assignment #1 ¹H 90° pw calibration

### Introduction

What is the 90 degree pulse width? The radio frequency pulse is described by its power and duration (time). Before the pulse the proton "spin" is oriented in the z direction and it cannot be observed in the spectrum. The rf pulse "flips" the spin into the xy plane so that it can be observed. The maximum signal is seen when the spin is completely in the xy plane without any z component. The length of time that this takes is called the 90 degree pulse. In the spectra below the length of the pulse is varied from 1 to 40 microseconds. The most accurate way to determine the 90 degree pulse width is to find the 360 and divide it by 4. This value is directly affected by the pulse power. As you increase the power, then the 90 degree pulse will decrease. Normally a power (tpwr) between 54-60 is used.



#### Instructions

1. Collect a ¹D 1H spectrum. Expand the region around one peak as shown in the spectrum below. Save this spectrum. File Name:



Take note of the previously calibrated 90° pw and power level and record them below. These values can be found on the acquisition page of the Acquire tab.

	Start Acquire	Process Show Tim	e Acquire	Stop	MoveSW Est. p	w90	Sequence	22	
	Default H1 Acquisition	PROTON		Display	Sequence	Arra	/5		
	Pulse Sequence	Data			Excitation				
	Channels	Spectral width	6410.3 Hz	<u>-</u>	Relaxation delay	1.000	sec 💌		
	Flags	Acquisition time	2.556 sec	-	First pulse	0.00	us 🔻	or O	degrees
	Future Actions	Complex points	16384		Inter-pulse delay	0.000	sec 💌		
3.	Overview	Scans Requested	8		Observe Pulse	3.30	us 🔻	or  45	degrees
		Completed			Receiver Gain:	30			_
ıt		Steady-State	0 ₪ ₪ off	$\leq$	Calibration: pw90	6.60	at Power	57	>
Ρ	Previously	calibrated	90°pw_		at po	wer_			<u>.</u>

2. Create an array of values for the pulse width (pw). Open the array pop-up window by clicking the arrays button found on the Default 1H tab of the Acquire panel. (Alternatively type array on the command line).

Start Acquire Default H1 Acquisition Pulse Sequence Channels Flags Future Actions	Process Acquire Acquire & Transform petral Width (ppm) 14 + -2 • Downfield 14.0 Upfield -2.0 Pulse Width (degrees) 45 • Enter pulse angle 45 Relaxation Delay (sec) 1 • Number of Scans 1 •	Show Time     Sequence     Arrays       Transform size:     32k     Image: Sequence       Line Broadening [Hz]     Image: Sequence     Sequence       Plotting turned off     Spectrum:     Image: Sequence       Parameters:     Image: Sequence     Image: Sequence       Image: Sequence     Image: Sequence     Image: Sequence	a >
- Temp 20.8 c	Spin Lock Probe Vite NutoX_DB_8790	Idle	ete

The array pop-up window is shown below:

		Arra	y Paramet	er		
Param Nam	e	Description		Size	Order	On/Off
A	Array Size:		Total	Time:0:0:1		
		UnArray		New A	rray	
Active Para	m:	None	Currer	nt Value:		
			1	Position		Value
Array Size:					INA	
First Value:						
Increment:						
Last Value:						
Inc. St	yle	None				
Rando	mize					
	Edit	Indo	0.6.0	Abandan		

Create the array by doing the following:

Type pw in the Param Name Box Click New Array Fill in the Array Size, First Value, and Increment Boxes Click enter on the keyboard.

The array needs to go past the  $360^{\circ}$  pw. This number can be estimated by taking the previously calibrated  $90^{\circ}$  pw value and multiplying it by 4. You will want your array to go past this point. In this example the  $90^{\circ}$  pw is 11.25 microseconds. The array must go past 11.25*4 (45 µseconds), therefore pw is arrayed from 1 to 55 µsecond in steps of 1 µsecond. After the array is created click the Close button at the bottom of the box.

		Arra	y Parameter	r		
Param Nam	ne	Description		Size	Order	On/Off
pw		Pulse width		55	1	On
	Array Size	:55	Total T	ime:0:2:54		
		UnArray		New A	rray	
Active Para	arn:	pw	Current	t Value:	5.95	
			1	Position	1 Va	lue
Array Size:		55	2		2 3	=
First Value		1	4		4	-
Increment:		1	6		6	
Last Value:		55	7		7	
Inc. S	tyle	Linear	9		9	
Rando	omize		10		10	
			12		12	
			13		13	
			14		14	

An appropriate delay, number of scans, absolute intensity must be now set and the experiment started. A delay of 2 seconds with 1 scan is chosen in the example. These are set by typing:

T--- N--- 20.00

d1=2 nt=1 ai ga on the command line.

<ul> <li>Applications Actions</li> </ul>		Tue Nov 20, 0
×	VnmrJ	
<u>F</u> ile <u>E</u> dit ⊻iew E <u>x</u> p	eriments Acquisition Process Tools Help	
: 🎽 📂 🔚 🚱 🛇	0	- 0 - L - M
		<b>3</b> 3×
exp2: Acquisition com	olete	
aph		
⊻ d1=2 nt=1 ai ga		
Exp:2 Seq: Proton	Index: 1	

When the experiment is completed display one spectrum from the array and phase it. In this example the third spectrum is displayed and phased. This is done by typing

ds(3) full aph on the command line.

	*	VnmrJ	=   =   ×
	<u>File</u> <u>E</u> dit <u>V</u> iew Experiment	ts <u>A</u> cquisition <u>P</u> rocess <u>T</u> ools <u>H</u> elp	
			<u></u>
$\bigcirc$	exp2: Experiment started exp2: Acquisition complete ds(3) full aph	$\sum$	
	Exp:2 Seq: Proton	Jirdex: 55	
		$\wedge \wedge \wedge \wedge$	
	2.92 2.90	288 2.86 2.84 2.82 2.80 2.78 2.76 2.74 2.77	2.70 2.68 ppm
	vp 12.0	cr vs delta 2.93 497.0	

To view the entire array and display the array values in the bottom text box, type dssh da on the command line. It may be necessary to adjust the vertical scale (vs) and the vertical position (vp) to make all spectra display completely on the screen. The text box with the array values is found in the Text Output tab of the Process panel.



The 90° pw is calculated by taking the 360° pw and dividing it by 4. The 360° pw is identified as the second null point in the curve. The numbers can be displayed

under each spectrum by typing dssl (not shown). In this example, the  $45^{th}$  spectra is the  $360^{\circ}$  pw. This corresponds to  $45 \,\mu$ seconds. Save the file. File name:



Print the spectra and array parameters by using the commands: pl('all') pssl page printon dg da printoff Label the printout with the 90, 180, 270 and 360 degree pulsewidth.

Record the  $90^{\circ}$  pw and power (tpwr): .

3. To find a more precise  $90^{\circ}$  pw create a second array that spans 4 µseconds around the  $360^{\circ}$  pw with 0.25 µseconds increments. In this example, the second array would be from 43 – 47 µseconds with an increment of .25 µseconds. Set up the array, collect the spectra and save the file.

Filenaname: ______.

Print the spectra and array parameters by using the commands: pl('all') pssl page printon dg da printoff

Record the 90° pw_____and power (tpwr):_____.

# A.4. 400 NMR Practical Assignment #2 gCOSY

Put your sample into the magnet. Lock, shim, tune and collect a proton experiment. Save the data and record the file name:

Calibrate a 90 degree pulse width. Record the calibrated value with the power at which is was collected: 90°pw:_____ tpwr:_____.

Create a new experiment by typing cexp(10) in the vnmrj command window. This will create experiment number 10. Go to experiment 10 by typing jexp10 in the command window.

Open the proton that was saved above. Set the 90 degree pulse width and the power to the calibrated value. You may fill in the boxes on the Acquisition page of the Acquire tab.

	Start Acquire	Process Show Tim	e Acquire	Stop	MoveSW Est. p	w90	Sequence	
<	<del>Default H1 Acquisition</del> Pulse Sequence Channels Flags Future Actions Overview	PROTON Data Spectral width Acquisition time Complex points	6410.3 Hz 2.556 sec 16384	Display	/ Sequence Excitation Relaxation delay First pulse Inter-pulse delay	Arra 1.000 0.00 0.000 0.000	ys sec ▼ us ▼ or 0 sec ▼	degrees
3. 		Scans Requested Completed Steady-State	8 0 0 🖌 off	<	-Receiver Gain: Calibration: pw90	30 6.60	at Power 57	

Minimize the sweepwidth by putting the cursors around the peaks and typing movesw. Collect a spectrum.

Record the sweepwidth (sw): _____.

Record the tuner offset frequency (tof): _____.

Save the file and record the file name: _____

Change the experiment into a gradient COSY experiment by using the Experiment drop down menu at the top of the page.



Click the gradient COSY button. It is found under Experiments>convert current experiments to do..> homonuclear correlation experiments>gradient COSY. Notice that the sequence at the top now reads gCOSY.



Look at the parameter list on the Overview page of the Acquire tab. Type dg to refresh this parameter list. Make sure that pw, tpwr, sw and tof are the same as you recorded above. If not, change them by typing the abbreviation and value (i.e. sw=4803). The sweepwidth for the second dimension should equal the sweepwidth for the first dimension in this experiment. Therefore sw1 should equal sw.

		A		i				8	i (
Start Acquire P	rocess	Show Time	Acqui	re Stop				Sequence	
 Defaults Acquisition Pulse Sequence PRESAT WET Channels Flags Future Actions Overview	ACQI seqfil sw at np ss d1 d2 nt ct 2D A( sw1 ni PRES/	UISITION gCOSY 6410.3 0.150 1924 32 1.000 0.020 1 0 0.020 1 0 0 0.020 1 0 0.020 1 0 0.020 1 0 0.020 1 0 0.020 1 0 0.020 1 0 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1 0.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.000 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020 1.020000000000	TRAN tn sfrq tof tpwr pw DEC dn dm SP temp spin gain	SMITTER H1 399.965 57 6.600 OUPLER C13 nnn ECIAL not used 0 30	GRAD gzlvlE gtE EDratio gstab hsglvl hsgt FL hs sspul SA date War solvent sample	IENTS 5102 0.001000 0.000500 6120 0.002000 AGS MPLE ch 8, 1993 d20	PRC sb fn 2D PF sb1 sbs1 proc1 fn1	CESSING -0.075 not used 2048 COESSING -0.020 not used 1p 2048	Text Basic Array Channels Solv.Supp Shims
	satmode wet Detai	e n n 1< in doss				II			Crear

Record the gain:

Type nt=1 ss=4 ni=1 ga

Watch the remote status box for a ADC overflow error. Raise the gain value until the signal overflows, then set the gain to 5 less than the overflow value. (The maximum value is 60).

Record the gain: _____.

Determine the necessary number of transients (nt) by observing how many transients are required to see the proton signal.

Record required nt: _____.

Type nt = required number as determined above, ss=32 ni=128 Type time to see how long the experiment will be.

.

Click the green Acquire button.

Save the experiment when it is completed. Record the filename:

# A.5. VNMRS-500 NMR testing information

*Equipment Components* Probes Triple Resonance Probe (TR_8064)



Varian High-Field Triple-Resonance Probe Manual, Publication: 01-999132-00E, page 7.



Varian High-Field Triple-Resonance Probe Manual, Publication: 01-999132-00E, page 5.



Figure 5. High-Field Triple-Resonance Probe Tuning Knobs

Varian High-Field Triple-Resonance Probe Manual, Publication: 01-999132-00E, page 8.

Indirect Detection Probe (id-8297)



Varian High-Field Indirect Detection Probe Manual, Publication: 01-999175-00E, page 6



Varian High-Field Indirect Detection Probe Manual, Publication: 01-999175-00E, page 4



Varian High-Field Indirect Detection Probe Manual, Publication: 01-999175-00E, page 8

Tune Interface Box

This box is located on the floor in front of the magnet. Set the channel to either 1, 2, or 3 as appropriate for the nuclei that is being tuned. Tune the probe until the readout is at least 50. The mtune program can also be used for tuning.



Remote Status Box:

The Remote Status Box sits on the table next to the computer. Do not leave the spin light flashing when you log off the computer. If the spin light is flashing, log back into VNMRJ and turn the spinner off.

XTC Ca are 🕘 xtg C XT VARIAN

FTS Sample Cooler

This is located behind the 500 console. It is use the sample temperature. The air flow must be 100 psi into compressor. Operation of this unit requires Var certification.



**Pneumatics Router** 

The router is located on the east wall of the NMR laboratory. Air flow is regulated here.



Gas flow sensor must be within green LED lights (adjusted within VNMRJ software)



Varian Pneumatics Router Manual, Publication: 01-999302-00B

# A.6. Transferring Data to your pc using Winscp

Data can be transferred using a sftp client. Winscp is freeware which can be downloaded from: <u>www.winscp.com/</u>

Host name: is the ip address of the computer which you would like to access. The 400 computer is 129.15.22.46. The User name and Password are your NMR username and password. See Section V of this manual for a list of IP addresses for the computers in the NMR facility.

Environment Directories SFTP SFTP SCP/Shell Connection Proxy Turnel SSH Authentication Bugs Ele protocol SFTP VIAItow SCP failback Bugs	Session Stored sessions Logging	Session Host name	Po <u>i</u> t number
SUP Software Ley file Prove Turnel SSH Key exchange Rule Software Ley file Ele protocol Bugs	Environment - Directories - SFTP	User name	Password
SH Key exchange Authentication Bugs	Connection Proxy Tunnel	Private <u>k</u> ey file	()
	SSH Key exchange Authentication Bugs	Protocol Eile protocol	SFTP V Allow SCP fallback
Select	Preferences		Select cold

After the host name, user name and password is entered, click Login to access the data. The left panel show the files on your computer, the right panel shows the files on the computer that you have logged into.

😼 vnmr1 - vnmr1@129	.15.12.236	- WinSCP						
Local Mark Files Comman	ids Session (	Options Remote H	lelp					
🚺 🛛 📦 • 🖙 📽	📀 🔤 🦧		∀ \$	Ø 🔿 Default	- 🚳 -			
*** == vowr1@129.15.15								
VIIII1@129.13.12		* <b>101</b>						
🛥 C: Local Disk 🛛 👻 🕓	[- + ]	🖭 🔝 🚮 💋	🔄 🗄	🗀 vnmr1	♥   (-+	- 🗈 🚺	🚮 😰	🗀 🗄
C:\Documents and Settings\su	isan\My Docur	nents		/home/vnmr1				
Name Ext ^	Size	Туре	Cha 🔨	Name 🔶 Ext			Size	Chang 📥
💼		Parent directory	5/6/	<b></b>				4/28/2
autotest_OU_3April2008		File Folder	4/3/	acrobat 🚞				8/24/2
AutoTripRes		File Folder	4/17	adobe 🚞				8/24/2
🚞 AutoTripRes_Murali		File Folder	4/17	🛅 .config				8/24/2
ChurchHistory		File Folder	3/3C	Carl eggcups				8/24/2
🚞 Joels Journal		File Folder	3/23	🚞 .gconf				5/7/20
📸 My Music		File Folder	4/15	🛅 .gconfd				5/7/20
My Pictures		File Folder	4/22	Canome .gnome				8/24/2
🕮 My Videos		File Folder	1/11	ignome2				4/28/2
i optimization assignment		File Folder	2/6/	gnome2_private				8/24/2
🚞 Personal_daniel_susan		File Folder	4/7/	.gstreamer-0.8				8/24/2
Receipts		File Folder	1/17	🛅 .kde				8/24/2
CRemote_operation_jan		File Folder	5/2/	Constant in the terminal interview of termin				8/24/2
methoxy_jan23_08_h1		File Folder	3/24	🛅 .mozila				8/25/2
🚞 steven 1. fid		File Folder	2/29	🛅 .nautilus				8/24/2
🖬 fj17_3	298,092	File	4/8/	Crhopenoffice1.1				12/11/:
🖬 methoxy_azo	612,901	File	3/24	🚞 .ssh				3/7/20
🔊 Thumbs.db	47,616	Data Base File	2/29	🛅 .thumbnails				8/29/2
Regular Grocery Store	162	Microsoft Word	1/19 🚩	🛅 . Trash				1/15/2
<			>	<				>
0 B of 8,449 KB in 0 of 90				0 B of 101 KB in 0 of 6	54			
/home/vnr	nr1\$							*
🖗 F2 Rename 📝 F4 Edit	B FS Copy	🚡 F6 Move 💣 F7	Create Di	rectory 🗙 F8 Delete	F9 Properties	👖 F10 Quit		
					A	SFTP-3	3	0:00:20
							uðu	

Find the file which you would like to transfer, highlight and click copy.

ୟ vnmr1 - vnmr1@129.	.15.12.236	- WinSCP				
Local Mark Files Comman	ds Session (	Options Remote H	-telp			
🌢 🛛 🖓 • 🕂 😫	🔉 📼 🎤		∀ ¢	🖉 🔿 Default 🔹 쨻 🗸		
			hadad halaad is			
vnmr1@129.15.12	•••••••••••••••••••••••••••••••••••••••	* <b>10</b>				
🖙 C: Local Disk 🛛 🗸		🖭 🔝 🚮 📝	😑 🗞	🔁 vnmr1 🛛 🗸 😓	- 🔝 🙆 🚮 🕼	😑 📴
C:\Documents and Settings\su	isan\My Docur	nents		/home/vnmr1		
Name Ext ^	Size	Туре	Cha 🔨	Name – Ext	Size	Chang 📥
<b>b</b>		Parent directory	5/6/	.openoffice-lock	0	12/11/:
autotest_OU_3April2008		File Folder	4/3/	a .recently-used	10,584	5/7/20
AutoTripRes		File Folder	4/17	.sversionrc	69	8/27/2
autoTripRes_Murali		File Folder	4/17	🖬 .viminfo	3,883	3/7/20
ChurchHistory		File Folder	3/30	.vnmrenv	1,552	9/7/20
🚞 Joels Journal		File Folder	3/23	.vnmrenv.bkup.070824.13:16	1,552	8/24/2
📸 My Music		File Folder	4/15	.vnmrenv.bkup.070907.16:39	1,552	9/7/20
📇 My Pictures		File Folder	4/22	👿 .vnmrenvsh	1,675	9/7/20
My Videos		File Folder	1/11	🖬 .vnmrenvsh.bkup.070824.13:16	1,675	8/24/2
optimization assignment		File Folder	2/6/	🖬 .vnmrenvsh.bkup.070907.16:39	1,675	8/24/2
🚞 Personal_daniel_susan		File Folder	4/7/	.vxresource	56	9/7/20
Receipts		File Folder	1/17	🖬 .vxresource.bkup	56	8/24/2
Remote_operation_jan		File Folder	5/2/	🖬 .vxresource.bkup.070824.13:16	56	8/24/2
methoxy_jan23_08_h1		File Folder	3/24	👿 .vxresource.bkup.070907.16:39	56	9/7/20
Consteven1.fid		File Folder	2/29	🖬 .Xauthority	696	5/7/20
🖬 fj17_3	298,092	File	4/8/	.xscreensaver	10,497	8/28/2
🖬 methoxy_azo	612,901	File	3/24	🖬 .zshrc	658	8/24/2
🔊 Thumbs.db	47,616	Data Base File	2/29	🔁 susan.ps	46,685	5/7/20
👼 ~\$qular Grocery Store	162	Microsoft Word	1/19 🚩			~
<			>	<		>
0 B of 8,449 KB in 0 of 90				0 B of 101 KB in 0 of 64		
/home/ym	nr1\$					*
P F2 Rename 7 F4 Evit	FS Copy	F6 Move 🎓 F2	7 Create Dir	ectory 🗙 F8 Delete 🐨 F9 Properties	F10 Ouit	
			2.0400 01			
				<u> </u>	SFTP-3 🗐	0:00:12

# A.9. Remote access to the 400 NMR via Putty and Real VNC

Putty and Real VNC are both freeware which can be downloaded from the internet.

#### **CONFIGURE PUTTY**

Open the Putty software

RuTTY Configuration	ı ? 🗙
Category:	
Session     Logging     Terminal     Keyboard     Bell	Basic options for your PuTTY session Specify the destination you want to connect to Host Name (or IP address) Port 22
Features     Features     Window     Appearance     Behaviour     Translation     Selection	Connection type: ○ <u>Raw</u> ○ <u>Lelnet</u> ○ <u>Rlogin</u> ○ <u>S</u> H ○ Serial Load, save or delete a stored session Savgd Sessions
Colours Connection Data Proxy Telnet Riogin BSH	Default Satings Load Vrimi Saye Delete
Serial	Close window on exit: Always Never Only on clean exit
ADOUT Help	<u>Upen</u> <u>Lancel</u>

#### Click on the SSH



#### Click on tunnel



#### Rutty Configuration **?**× Category: Options controlling SSH port forwarding 🚊 Terminal ^ - Keyboard Bell Features Port forwarding Local ports accept connections from other hosts Remote ports do the same (SSH-2 only) - Window Appearance Behaviour Forwarded ports: <u>R</u>emove Translation - Selection - Colours Add new forwarded port: Connection Data Proxy Telnet Source port 5901 Add Destination localhost:5901 O Remote Rlogin ⊡ SSH Local O Dynami 💽 Aut Kex Auth TTY X11 Tunnels ~ Bugs <u>O</u>pen <u>C</u>ancel <u>A</u>bout <u>H</u>elp

#### Put 5901 into the Source port and localhost:5901 into the destination

#### Click Add



Click Session (it is at the very top)


Put the IP address of the computer that you wish to connect to into Host Name and a name in the saved Session and click Save. The IP address of the 400 NMR is 129.15.22.46. In this case, I saved it as "400 NMR".



#### Using PUTTY:

Open the Putty software. Highlight your saved Session and Click Load.



Click Open



#### A terminal window will appear



Log in with your NMR login and password

📴 129.15.12.236 - PuTTY	
login as: rhc	<u>^</u>
Incel29.13.12.230 B password:	

Just return when prompted for the input display server name.



Type vncserver at the prompt. The session number must be 400Computer:1 or the vncviewer won't open in the next step.

400Computer~vnmr1 :/home/vnmr1	
al	
login as: vnmr1	
e vnmr1@129.15.22.46's password:	
Last login: Tue Jan 25 14:43:34 2011 from d-ip-10-196-225-86.wi	reless-rc-priv.no
r.ou.edu	
xrdb: Can't open display ''	
vnmr1 1>vncserver	
New '400Computer:1 (vnmr1)' desktop is 400Computer:1	
Starting applications specified in /home/vnmrl/.vnc/xstartup	
Log file is /home/vnmr1/.vnc/400Computer:1.log	
vnmr1 2>	
	· · · · · · · · · · · · · · · · · · ·

Open vnc viewer

type localhost:5901 in the Server Blank. Then click OK

VNC View	er : Conne	ction Details	
VO	Server:	localhost:5901	~
<u>VC</u>	Encryption:	Always Off	~
About	. Optic	ons OK	Cancel

Optional: If you want to share the screen, then you must click Shared connection under the Misc tab of the Options menu



Put your remote password in the password blank

1	VNC View	er : Auther	ntication [No Encryption]	
ļ		Username:		ОК
		Password:		Cancel

This will open the 400 desktop.

Run your experiments.

When you are finished, log out of the desktop under actions at then top. Then type pkill vncserver pkill vnc at the command prompt. This is extremely important. These are the commands that will reset the session number to :1.



# A.8. Glossary of Common NMR Commands and Terms

aa	abort acquisition, hard stop
acqi	Open acqi window if the button has been lost
ai	absolute intensity mode
alock	autolocking routine, <b>alock='y'</b> for autolocking, <b>alock='n'</b> for typical manual locking
aph	autophasing, not recommended for most spectra
array	macro for setting up an arrayed experiment
<b>at</b> points ( <b>np</b> )	acquisition time, set by spectral width ( <b>sw</b> ) and number of data
axis='p'	specify ppm or hertz for the axis, i.e. <b>axis='p'</b> or <b>axis='h',</b> use
axis= pu	referencing in an indirect detect experiment with >decref
BPsvf	BioPack save file command: saves all pieces associated with an Experiment including any shape pulses, global file, probe file, etc. Usage: >BPsvf('filename')
bc	baseline correct
bs	block size, data is stored to the disk every time an increment of bs is reached, i.e. <b>bs=16</b> , every 16 scans data is saved and can be transformed
<b>cd</b> directory	change directory, changes the directory back to home default
cexp(#)	create the experiment #
COSY	correlation spectroscopy, a 2-D experiment, homonuclear one-bond J coupling
center	re-size a 2-D spectrum to a centered square, same as [DispMenu] [Size][Center]
ct	completed transients
CZ	clear all integral reset points

**d1** delay time between scans, required to allow for T1 relaxation, in seconds

da	display arrays
dc	drift correct
dconi	display interactive color map (2-D)
df	display FID
dfrq	decoupler frequency (2 nd channel)
dg	display first text screen of parameters
dli	display list of integrals
dm	decoupler mode
dmf	decoupler modulation frequency
dmm	decoupler modulation mode (c, g or w)
dn	decoupler nucleus (2 nd channel)
dof	decoupler offset (Channel 2 transmitter offset) in Hertz
dp	double precision, set to 'y'
dpcon	display contours in 2-D spectra
dpf	display peak frequencies
dpir	display integrals on screen (requires <b>vp=12</b> )
dps	display pulse sequence
dpwr	decoupler power
dres	display digitial resolution
ds	display spectrum
dscale	display the scale (in ppm or Hertz)

dsn	display signal to noise
dss	display stacked spectra
dssa	display stacked spectra automatically
dssh	display a series of spectra in an arrayed experiment
dssl	display corresponding numbers in arrayed experiment
explib	display experiment library (or list of current exps.)
f	display the whole spectrum
fb	filter bandwidth
fn/fn1	Fourier number for direct (fn) detected dimension, indirect (fn1)
foldt	fold COSY type spectrum along diagonal
ft	fourier transform the data (no weighting functions)
full	display over the whole screen
fullt	re-size a 2-D spectrum for full with traces, same as [DispMenu][Size][Full with Traces]
ga	get acquisition (start acquisition and transform data)
gain	the receiver gain, to see the value of <b>&gt;gain?</b> to set the value <b>&gt;gain=40</b>
gCOSY	gradient COSY, same as a COSY, fewer scans required
gettext	will bring up a very simple editing window to type text into
gf	gaussian weighting function
gmapsys	start gradient shimming routine and open menu
go	acquire the spectrum, don't transform
gzsize	number of Z shims to use in gradient shimming
НМВС	heteronuclear multiple bond coherence, 2-D heteronuclear experiment

**HMQC** heteronuclear multiple quantum coherence, 2-D heteronuclear experiment,

One-bond correlations

- jexp# join a particular experiment, jexp2
- **Ib** line broadening weighting function (exponential)
- left set display to left side of screen
- **LOCK** The deuterium nuclei in the sample are used to maintain a "lock" on the sample. The nuclei are used to monitor and correct for any drift in the magnetic field. If the field "drifts" or changes in strength, the precessional frequency of a nucleus will change accordingly. In a pulsed lock system, the field is monitored by observing the resonance frequency of the deuterium nucleus of the solvent (i.e. D2O). The resonance frequency of the nucleus is compared to a reference frequency in the spectrometer and any changes are corrected by adjusting Z0.
- **Lock gain** the amplification of the deuterium NMR signal, increases the size of the signal, but also increases any other signals or noise that may be present.
- **Lock phase** the phase angle used to control the phase of the deuterium NMR signal and the phase of the reference signal for the deuterium lock, normally needs very little if any adjustment.
- Lock power The quantity of rf energy used to irradiate the deuterium nucleus, controls the amplitude of the rf pulse at deuterium frequencies. Must be large enough to produce a signal for the deuterium but still below the saturation limit. If the power is too high, the lock signal may decrease in intensity.
- **Ip** left phase, first order phase correction
- **man** a very useful command to access the manual on an experiment i.e.>man('noesy')
- **movesw** move sweep width, first enclose the region for the sweep width with the cursors, then type **movesw**, this will move the **tof**
- **movetof** move transmitter offset, place cursor on peak or position to the center point of the spectrum, type **movetof**, does not change **sw** value.

ni	number of increments
nl	nearest line
np	number of data points acquired in the FID
nt	number of transients or scans
p1	another pulse that can be used in certain experiments
pad	pre-acquisition delay
page	sends plotting commands to the printer
рар	print parameters on plot, long version
phase	used to set phase selection in multi-dimensional experiments
phase(180)	phase the spectrum - 180 degree flip
pl	plot the spectrum
plfid	plot the FID
pll	plot line list
pltext	plot the text
plww	plot arrayed spectra in whitewash mode
рра	print parameters, written out on plot
ppf	print peak frequencies
printon/ printoff	starts the printer job and ends the job
pw	pulse width measured in microseconds
pw90	the 90-degree pulse width, corresponds to the amount of time the transmitter is on in order to achieve a 90 degree tip angle
ra	resume acquisition stopped with <b>sa</b>
rl(4.6p)	reference a line to 4.6 ppm

rl1(4.6p)	reference a peak in a 2-D homonuclear experiment in f1 to 4.6 ppm
rl1(77d)	reference a peak in 2-D indirect detection in f1 to 77 ppm
rp	right phase, zero-order phase correction
rts	retrieve shims
sa	stop acquisition, this is a soft stop which means it will stop after the next FID
sb	sinebell weighting function
sd	set decoupler
sfrq	spectrometer frequency in MHz
SHIM field	The process of "shimming" a sample is to minimize or eliminate any
neid	differences across a sample. Eliminating these differences will lead to narrower lines and increased intensity.
SS	steady state scans, scans put in before acquisition really begins to create a steady state
su	set up the experiment, must be used when retrieving shims, setting nucleus for tuning, changing the temperature, etc.
svf	save file
svp	save parameters
SVS	save shims
SW	the spectral width used to sample NMR signals, directly related to the chemical shift range for a given nucleus, given in hertz, sets the rate at which data is sampled.
tn	transmitter nucleus (i.e. H1, channel 1 on tune box)
TOCSY	total correlation spectroscopy, 2-D homonuclear proton experiment, through bond couplings, multiple bonds
tof	transmitter offset (Channel 1) in Hertz

tpwr	transmitter power in dB
trace	mode for 2-D or greater display (trace='f1' or trace='f2')
TUNE	Tuning a sample reduces the amount of power reflected back to the transmitter
vp	vertical position
vs	vertical scale
vs2d	vertical scale for a 2-D spectrum
vsadj	vertical scale adjust, adjusts to tallest peak in display
<b>vttype</b> changes,	setting for temperature control, vttype=2 allows temperature
	vttype=0 does not allow temperature changes
wbs('wft')	with the next block store, transform the data
wft	weighted Fourier transform
wti	open interactive weighting
Z0	The Z0 allows the operator to match the resonance frequency of the deuterium to the reference frequency for the deuterium lock.

## A.9. A few useful linux commands

To open a terminal window on facility workstations: Hold down the right mouse button on the wallpaper of the screen. The menu will appear. Select "Programs" and then select "Terminal".

Commands that can be used in the terminal window: Note: .fid files are directories in UNIX

## cd - change directory

>cd
goes back to the home directory
>cd ..
go up one directory to the parent directory
>cd /net/inova6001
for changing directories from one facility workstation to
another facility workstation

### compress - compress a file to take up less disk space

#### >compress data.tar

make a tar file (see tar command) from a .fid directory and then compress the data using this command. This results in a file like: data.tar.Z

to reverse the compression:

### >uncompress data.tar.Z

**cp** - copy files

>cp thisfile thatfile
copies the contents of "thisfile" into a file called "thatfile"
>cp -r directory1 directory2
same thing only for directories
>cp filename directory1/subdirectory1/.
copy a file into another directory's subdirectory

du - disk usage

>du -k mydirectory

Be in the parent directory of the directory to be checked. Put in the username (usually the user's main directory) and can check the usage on the whole account. Otherwise, also put in a subdirectory name and check the amount of disk space it takes up.

exit - exits the terminal window >exit Is - list files and directories

>ls

lists all of the files and subdirectories in the current directory >Is -I

lists all of the files and subdirectories with dates and other info.

## mkdir - make directory

### >mkdir newdirectory

creates a new directory called "newdirectory"

**more** - used to displaythe contents of a file without editing mode >more myfile

displays the text of the file, page ahead with the space bar

- mv move files
  - >mv file1 directory2/.

moves file1 into directory2 and calls it the same name(indicated by the period)

## passwd - change the account password

### >passwd

Terminal will prompt for the old password and the new password (twice). Use different cases and numbers in the password.

## **pwd** - print working directory

## >pwd

shows the current directory location

## rm - remove files

>rm junkfile
removes a file to be deleted
>rm -r junkdirectory
>rm -r junk.fid
removes a directory

rmdir - remove directory
>rmdir directory3
>rmdir junk.fid
.fid files are actually directories and have to be removed this way

tar - tape archive - This command allows packing a directory into a "tar" file which acts like a single file. So, data can be stored as one directory, pack it into a tar file and move it around like a single file. Then unpack it into the new location >tar -cvf directoryname.tar directoryname

>tar -cvf data.fid.tar data.fid
Creates a .tar file from the directory, "directoryname"
>tar -xvf directoryname.tar directoryname
>tar -xvf data.fid.tar data.fid
Extracts the .tar file back into a normal directory.

**tcsh** - changes the type of terminal to one that will allow actions like up arrow for the last command, an enhanced c-shell.